

# THINK POSSIBLE

## **BioTek Microplate Instrumentation**




Think Possible


BioTek


BioTek Instruments, Inc., headquartered in Winooski, VT, USA, is a worldwide leader in the design, manufacture, and sale of microplate instrumentation and software. For over 45 years, our products have been designed and manufactured in the U.S.A. BioTek's instrumentation is used to aid in the advancement of life science research, facilitate the drug discovery process, provide rapid and cost-effective analysis and to enable sensitive and accurate quantification of a wide range of molecules across diverse applications.

Our company-wide commitment to quality and value is backed by superior customer care, technical service centers, scientific application experts, and a knowledgeable sales force. Our commitment and focus helps to ensure your processes will be rapid, efficient and successful.

This catalog provides an overview of our complete line of microplate instrumentation, and features Cytation™, the first instrument to combine automated digital microscopy and conventional microplate detection. For more detailed information and product specifications, visit our web site at www.biotek.com.







NH<sub>2</sub>

| ▶ Imaging & Microscopy                            |          |
|---|----------|
| Cytation™   |          |
| Gen5™ Data Analysis Software                      | 8        |
| ▶ Detection                                       | 10       |
| Hybrid Technology Multi-Mode Microplate Readers   | 4.0      |
| , 5,  |          |
| Synergy H1  | 14       |
| Multi-Mode Microplate Readers                     |          |
| Synergy 2   |          |
| Synergy HTX                                       |          |
| Multi-Mode Reader Comparison Chart                | 20       |
| Fluorescence Microplate Readers                   |          |
| FLx800™   | 22       |
| Absorbance Microplate Readers                     |          |
| PowerWave <sup>™</sup> HT                         | 23       |
| Epoch™ 2  | 24       |
| Epoch   | 26       |
| ELx808™   | 28       |
| ELx800 <sup>™</sup>                               | 29       |
| Absorbance Reader Comparison Chart                |          |
| Gen5 Data Analysis Software                       |          |
| Take3™ Micro-Volume Plate                         |          |
| Reader Accessories                                |          |
| N 122-4 11  | ٠,       |
| Liquid Handling                                   | 30       |
| Combination Microplate Washer Dispenser<br>EL406™ | 20       |
|   | 30       |
| Microplate Washers                                | 40       |
| 405™ Touch  |          |
| ELx405™ Select Deep Well                          |          |
| ELx50™  |          |
| Washer Comparison Chart                           | 44       |
| Microplate Dispensers                             |          |
| MultiFlo™ FX                                      |          |
| MicroFill™  |          |
| Dispenser Comparison Chart                        |          |
| Liquid Handling Control Software                  |          |
| Washer & Dispenser Accessories                    | 53       |
| Microplate Pipetting Systems                      |          |
| Precision <sup>™</sup>                            | 54       |
| ▶ Robotics  | 56       |
| NEW BioSpa™ 8 Automated Incubator                 |          |
| BioStack™ Microplate Stacker                      |          |
|   |          |
| ► Service & Support                               | 62       |
| Compliance  | 64       |
| Applications Support                              | 66       |
| Global Service and Support                        |          |
| • •   |          |
| Decades of Innovation                             | 68<br>71 |
| COMACT IMORNATION                                 | / !      |





"It is very easy to use & produces high quality images. Useful for both IF of tissues and cells from 96 well plate to petri dishes."

(Cytation 3 Cell Imaging Multi-Mode Reader)



The fusion of automated digital widefield fluorescent microscopy and multi-mode microplate detection in the Cytation Cell Imaging Multi-Mode Readers offers life science researchers a unique and powerful tool. Whether studying fixed or live cells, tissues or whole organisms, Cytation's automated fluorescence microscopy and conventional multi-mode detection of a wide array of fluorescent probes provides quantitative phenotypic data from images. These cellular images provide rich visualization and multi-parametric information concerning many biological processes.

This unique, patent pending design results in workflow efficiency and a reduction in data analysis and storage requirements, making Cytation the ideal instrument as a standalone automated imager or in combination with microplate reading as a total system.

### Cytation™ Cell Imaging Multi-Mode Reader

Cytation™ is a uniquely integrated, configurable system that combines automated digital widefield microscopy with conventional multi-mode microplate detection to provide phenotypic cellular information and well-based quantitative data. With up to 60x magnification, the microscopy module provides high-quality cellular and subcellular visualization. The multimode detection module provides high quality quantitative and qualitative data in all detection modes. All controlled with Gen5<sup>™</sup> software, specifically designed for uncomplicated processing of even the most complex assays

## Powerful Automated Digital Microscopy

The Cytation family includes Cytation 5 and Cytation 3. Each is available in multiple upgradable configurations; from basic microscopy to complex image collection and analysis in fluorescence, brightfield, color brightfield and phase contrast, Cytation offers critical methods like single- and multi-color, time lapse, montage and z-stacking. Available laser autofocus and image-based autofocus ensure fast and accurate image aquisition with minimal phototoxicity across a broad application range. Gen5 software offers unique features like user-trained autofocus and available joystick control for truly automated imaging with the highest quality results. Powerful image processing like image stitching, z-projection and digital phase contrast, are all available in this cost effective system.

#### Live Cell Assays

To create the ideal environment for live cell assays, Cytation includes 4-Zone<sup>™</sup> incubation up to 65 °C and a gas controller to monitor and control CO<sub>2</sub> and O<sub>2</sub> levels in the system. Linear, orbital

and double orbital shaking help keep cells gently agitated or well suspended to optimize many cell based assay protocols. To fully automate live cell assay workflows, Cytation integrates with BioSpa™ 8 Automated Incubator. From sample prep to image analysis, BioSpa 8 offers walkaway automation.

#### **Hybrid Technology**

BioTek's patented Hybrid
Technology, available with the
multi-mode detection modules
for Cytation, combines high
performance filters with variable
bandwidth monochromators,
providing convenience, versatility
and excellent performance.
Luminescence, UV-Vis absorbance,
time-resolved fluorescence,
fluorescence polarization and
Alpha detection modes greatly
increase the application range of
the system.

- ▶ 2D and 3D cell imaging and analysis
- ► Cell proliferation studies
- Cytotoxicity
- Biomarker quantification
- Drug discovery
- Genetic analysis
- Drug absorption and metabolism
- ▶ Biologics drug discovery and development
- ► Environmental testing
- Food safety
- Nucleic acid quantification
- Protein quantification







### **Specifications**

| General                           | Cytation 5  | Cytation 3   |  |
|-----------------------------------|---|--|--|
| Detection modes                   | UV-Vis absorbance<br>Fluorescence intensity<br>Luminescence<br>Fluorescence polarization<br>Time-resolved fluorescence<br>Alpha                     | UV-Vis absorbance<br>Fluorescence intensity<br>Luminescence<br>Fluorescence polarization<br>Time-resolved fluorescence                         |  |
| Read methods                      | Endpoint, kinetic, spectral scannir   | ng, well area scanning   |  |
| Microplate types                  | Monochromator: 6- to 384-well plates<br>Filters: 6- to 1536-well plates<br>Imaging: 6- to 1536-well plates  |  |  |
| Other labware supported           | Microscope slides, Petri and cell<br>(T25), counting chambers (hemor<br>Take3 Micro-Volume Plates   | Microscope slides, Petri and cell culture dishes, cell culture flasks<br>(TZ5), counting chambers (hemocytometer)<br>Take3 Micro-Volume Plates |  |
| Temperature control               | 4-Zone™ incubation to 65 °C with Condensation Control™  | 4-Zone incubation to 45 °C with Condensation Control™  |  |
| Shaking                           | Linear, orbital, double orbital   |  |  |
| Software                          | Gen5™ Data Analysis Software<br>Gen5 Secure for 21 CFR Part 11  | compliance (option)  |  |
| Automation                        | BioStack <sup>™</sup> and 3rd party automa<br>BioSpa <sup>™</sup> 8 Automated Incubator   | tion compatible<br>compatible  |  |
| Gas Control                       | Range: 0 - 20% (CO <sub>2</sub> ); 1 - 19% (C   | D <sub>2</sub> )   |  |
| Light Source                      | Xenon flash: FL and Abs<br>100 mW laser: Alpha  | Xenon flash: FL and Abs  |  |
| Imaging System                    |   |  |  |
| Imaging mode                      | Fluorescence, Brightfield<br>Color brightfield<br>Phase contrast  | Fluorescence<br>Brightfield  |  |
| Imaging method                    | Single color, multi-color, montag   | e, time lapse, z-stacking  |  |
| Image processing                  | Z-projection, digital phase contrast, stitching   |  |  |
| Camera                            | 16-bit gray scale, Sony CCD, 1.2<br>0.3 µm/pixel at 20 x  | 5 megapixel.   |  |
| Objective capacity                | 6 user-replaceable objectives   | 2 user-replaceable objectives  |  |
| Objectives available              | 1.25x, 2.5x, 4x, 10x, 20x, 40x, 60  | )x   |  |
| Phase objectives available        | 4x, 10x, 20x, 40x   |  |  |
| Imaging<br>cubes/capacity         | 4 onboard, user-replaceable cub   | es; >15 colors available   |  |
| Imaging LED cubes                 | Available: 365 nm, 390 nm, 465 623 nm, 655 nm, 740 nm   | nm, 505 nm, 523 nm, 590 nm,  |  |
| Automated functions               | User-trained autofocus, autofocus, autoexposure, auto-LED intensity   |  |  |
| Autofocus method                  | Image-based autofocus; Laser au   | itofocus (option)  |  |
| Positional controls               | Software control<br>Joystick control (option)   | Software control   |  |
| Image collection rate             | 96 wells, 1 color (DAPI), 4x: 6 m<br>96 wells, 3 colors, 4x: 12 minutes   | inutes (3 minutes laser AF)<br>s (7:30 minutes laser AF)   |  |
| Image analysis<br>software option | Gen5 Image+: Advanced image<br>Gen5 Secure Image+: Adds 21 C  | analysis<br>CFR Part 11 features   |  |
| Absorbance                        |   |  |  |
| Wavelength selection              | Monochromator   |  |  |
| Monochromator                     | Bandwidth: 4 nm (230-285); 8 nm<br>Wavelength accuracy: ± 2 nm<br>Repeatability: ± 0.2 nm<br>Range: 230 to 999 nnm, in 1 n<br>Resolution: 0.0001 OD |  |  |
|                                   |   |  |  |

|                           | Cytation 5  | Cytation 3                   |
|---------------------------|---|------------------------------|
| Optical density           | Accuracy: <1% at 2.0 OD; <3%<br>Linearity: <1% from 0 to 3.0 OI<br>Repeatability: <0.5% at 2.0 OI<br>Stray light: 0.03% at 230 nm | )                            |
| Pathlength correction     | Yes   |                              |
| Reading speed             | 96 wells: 11 s, 384 wells: 22 s (k  | kinetic interval)            |
| Fluorescence Intensity    | •   |                              |
| Wavelength selection      | Quad monochromators (top/bo<br>Filters (top)  | vttom)                       |
| Wavelength range          | Monochromators: 250 - 700 nm<br>Filters: 200 - 700 nm (850 nm o   | n (850 nm option)<br>ption)  |
| Monochromator bandwidth   | Variable; from 9 to 50 nm, in 1 nm increments   | Fixed, 16 nm                 |
| Dynamic range             | 7 decades   | 5 decades                    |
| Sensitivity (Fluorescein) | Filters: 0.25 pM (0.025 fmol/well<br>Quad Monochromator: 2.5 pM (<br>384-well plate) - top<br>4 pM (0.4 fmol/well, 384-well pla   | 0.25 fmol/well,              |
| Reading speed             | 96 wells: 11 s, 384 wells: 22 s (k  | kinetic interval)            |
| Luminescence              |   |                              |
| Wavelength range          | 300 - 700 nm  |                              |
| Dynamic range             | > 6 decades   |                              |
| Sensitivity               | Monos: 20 amol ATP (flash)<br>Filters: 10 amol ATP (flash), 100   | amol (glow)                  |
| Fluorescence Polarizat    | tion  |                              |
| Wavelength selection      | Filters   |                              |
| Wavelength range          | 280 - 700 nm (850 nm option)  |                              |
| Sensitivity               | 1.2 mP standard deviation at 1  | nm fluorescein               |
| Time Resolved Fluores     | scence  |                              |
| Detector                  | PMT   |                              |
| Wavelength selection      | Quad monochromators (second Filters (top)   | dary mode)                   |
| Wavelength range          | Filters: 200 - 700 nm (850 nm o   | ption)                       |
| Sensitivity               | Filters: Europium 40 fM (4 amol/well, 384-well plate)<br>Monos: Europium 1200 fM (120 amol/well, 384-well plate)                  |                              |
| Alpha Detection           |   |                              |
| Light source              | 100 mW 680 nm laser   |                              |
| Wavelength selection      | Filters (top)   |                              |
| Sensitivity               | 100 amol bio-LCK-P<br>(384-well low volume plate)   |                              |
| Reagent Dispensers        |   |                              |
| Number                    | 2 syringe pumps   |                              |
| Supported labware         | 6- to 384-well plates, Petri dish   | es                           |
| Dead volume               | 1.1 mL, 100 μL with back flush  |                              |
| Dispense volume           | 5-1000 μL in 1 nm increments  |                              |
| Physical Characteristic   | :s  |                              |
| Power                     | 250 Watts max.  |                              |
| Dimensions                | 20" D x 16.5" W x 17.5" H (50.  | 8 cm x 41.91 cm x 44.5 cm)   |
| Weight                    | 80 lbs (36.3 Kg)  |                              |
| Regulatory                |   |                              |
| Regulatory                | CE and TUV marked. RoHS Cor<br>Diagnostic use are available.  | mpliant. Models for In Vitro |

### Gen5™ Data Analysis Software for Cytation™

Gen5™ automates the entire imaging process: all the steps of image acquisition are fully automated including X, Y and Z movements, focus, LED intensity and camera exposure. Image analyses like cell counting and subpopulation analysis happen on the fly as images are acquired; Gen5's data analysis tools: kinetic analysis, curve fitting, EC50, potency calculation and others can then be applied to image analysis results. With Gen5 you can go from samples to final results in one click.

#### **Imaging Made Easy**

Gen5 was designed for users without prior automated microscopy experience. Image capture and acquisition is done through a simple interface that includes auto focus and auto exposure algorithms. Image

analysis can be as simple as a one click process with visual feedback to fine-tune analysis parameters.

#### **Automatic Cell Counting**

Gen5 can automatically count cell nuclei in samples. In a 96- or 384-well microplate, this means performing a count on tens of thousands of cells in a matter of minutes. This makes Cytation a very powerful tool for cellular assays where multiple conditions need to be tested and changes in cell population need to be monitored.

#### **Sub-Population Analysis**

Gen5 can sort cells by features such as intensity, or morphology (size, perimeter, circularity). This enables applications such as transfection efficiency, nuclear translocation or cell cycle assays where multiple cell sub-populations are present in the samples.

#### Hit Picking

Gen5 can be programmed to quickly scan the plate with the standard plate reader optics and hit-pick the wells that meet a defined threshold for imaging. This saves valuable time, reduces data collection, analysis and storage requirements – ultimately reducing costs.

#### Advanced Image Processing

Cytation Cell Imaging Readers can collect a montage of images from live cell kinetics, tissue sections and other large objects or samples. Gen5's powerful stitching capability seamlessly creates a complete picture of the sample with great accuracy, providing more meaningful data. After acquisition, Gen5 offers multiple processing tools prior to image analysis; automatic stitching of montages, z-projection of z-stacks, and digital enhancement of brightfield images.

#### **Key Features:**

Powerful instrument control


- Reliable autofocus methods including unique user-trained autofocus, and laser autofocus option
- Automatic camera gain, exposure and LED intensity settings
- ► End-point, montage, Z-stack and time-lapse read modes

#### Image pre-processing tools

- ▶ Built-in hot-pixel correction
- Automated image pre-processing (flattening, smoothing, background correction)
- ▶ Image stitching, Z-projection
- Digital phase contrast algorithm

#### Image and data analysis tools

- ▶ Automatic cell-counting and confluence
- Powerful sub-population analysis
- Image statistics
- ► Full data analysis software (EC50, standard curves, kinetic analysis and more)




HCT 116 tumoroids at 20x. Z-stacked image




Whole mouse at 2.5x. Stitched montage, color brightfield









Zebrafish at 4x. Texas Red (blood cells) & GFP (vasculature)




BPAE cells at 20x. Texas Red, Phalloidin & GFP



Euglena gracilis at 60x. Brightfield & chlorophyll (Red)



Lillium anther at 20x. Phase contrast, GFP & Texas Red




HT-1080 cells, migration assay at 2.5x, 1 hr to 12 hr run





"Epoch is a nice plate reader that can read microplates of various formats, up to 384 wells. It is a pretty robust instrument, we've had it for three years now with frequent use, hasn't needed any service. The instrument can do wavelength scans as well. Software (Gen5) is easy to use, and also very versatile... I'd recommend this instrument for any lab that needs a UV-Vis microplate reader."

(Epoch Microplate Spectrophotometer)



BioTek offers an extensive range of microplate readers, from the Synergy™ Neo2 Multi-Mode Microplate Reader to the ELx800™, a basic ELISA reader used in tens of thousands of laboratories around the globe. Included in the BioTek reader product range are Hybrid readers, multi-mode readers, fluorometers, luminometers and a variety of both monochromator-based spectrophotometers and filter-based absorbance readers.

For assays like nucleic acid and protein quantification, where very small sample size is critical, BioTek's Take3 $^{\text{TM}}$  Micro-Volume Plate offers the ability to measure multiple samples as small as 2  $\mu$ L in monochromator-based multi-mode and absorbance readers. BioTek's microplate readers come with the powerful Gen5 $^{\text{TM}}$  Data Analysis Software and are compatible with BioStack $^{\text{TM}}$  and many third-party automation products, to provide increased throughput and unattended operation. To automate live cell assays, several readers are compatible with the new BioSpa $^{\text{TM}}$  8 Automated Incubator.

### Synergy™ Neo2 Multi-Mode Reader

Synergy™ Neo2 Multi-Mode
Microplate Reader is the most
advanced, high performance
multi-mode microplate reader
available today. The outstanding
features of Synergy Neo2 allow
complex applications including
both biochemical and cellbased assays, to be performed
rapidly, efficiently and with
uncompromised performance, in
all detection modes.

#### Patented Hybrid Technology™

Some workflows benefit from the flexibility of monochromator-based optical systems; there's no need to purchase multiple filters, and when a fluorophore's spectral peaks are unknown, monochromators can scan to find the ideal excitation and emission peaks. Other assays require the high sensitivity found in filter-based optical systems. BioTek's patented Hybrid Technology offers both major

benefits in a single platform, so there's no compromise of performance or flexibility.

### Variable bandwidth quad monochromators

Synergy Neo2's monochromators, have variable bandwidths for excitation and emission. Selectable from 3 nm to 50 nm in 1 nm increments, these continuously variable bandwidths help optimize detection of some fluorophores. Detection parameters for complex multi-plexed assays like FRET and SNPs can be fine-tuned for the highest signal with the lowest crosstalk – and the results you expect.

#### Ultra-fast plate processing speeds with multiple PMT detectors

High throughput isn't just about fast plate reading – a high

throughput multi-mode reader should handle common and complex assays with equally high performance even in 1536-well plates. Synergy Neo2 has dual PMTs for top measurement, so FP, FRET, TR-FRET and other ratiometric measurements are processed quickly and with excellent results. Up to four PMTs are available in Synergy Neo2, for the greatest speed and flexibility.

## Controlled environment for live cell assays

Along with incubation to 65 °C and shaking, Synergy Neo2 can be equipped with a CO₂/O₂ controller to provide the ideal environment for robust live cell assays. Direct bottom detection provides ultra sensitivity for measuring cell-based fluorescence intensity. To automate live cell workflows, Synergy Neo2 integrates with BioSpa™ 8 Automated Incubator.

- ▶ HTS screening
- Drug absorption and metabolism
- Biologics drug discovery and development
- Drug discovery
- Cell proliferation
- Cytotoxicity
- ▶ Biomarker quantification
- Genetic analysis
- Environmental testing
- Food safety
- Nucleic acid quantification
- Protein quantification



Synergy Neo2 shown with optional high-speed microplate stacker





### Specifications:

| General                                    |   |
|--|---|
| Detection mode                             | UV-Vis absorbance<br>Fluorescence intensity<br>Luminescence<br>Fluorescence polarization,<br>Time-resolved fluorescence<br>Alpha      |
| Read methods                               | Endpoint, kinetic, spectral scanning, well area scanning  |
| Microplate types                           | 6- to 1536-well plates  |
| Other labware                              | Petri and cell culture dishes<br>Take3™ Micro-Volume Plates   |
| Temperature control                        | 4-Zone <sup>™</sup> incubation to 65 °C with Condensation Control <sup>™</sup>  |
| Shaking                                    | Linear, orbital, double orbital   |
| Software                                   | Gen5 <sup>™</sup> Data Analysis Software<br>Gen5 Secure for 21 CFR Part 11 compliance (option)  |
| Automation                                 | BioStack™ and 3rd party automation compatible<br>BioSpa™ 8 Automated Incubator compatible   |
| CO <sub>2</sub> and O <sub>2</sub> control | Range: 0 - 20% (CO <sub>2</sub> ); 1 - 19% (O <sub>2</sub> )  |
| Barcode reader                             | 1D and 2D camera-based scanner  |
| Absorbance                                 |   |
| Light source                               | Xenon flash   |
| Detector                                   | Photodiode  |
| Wavelength selection                       | Monochromator   |
| Wavelength range                           | 230 to 999 nm, in 1 nm increments   |
| Monochromator bandwidth                    | 2 nm (230-285 nm); 4 nm (>285 nm)   |
| Dynamic range                              | 0 - 4.0 OD  |
| Resolution                                 | 0.0001 OD   |
| Pathlength correction                      | Yes   |
| Monochromator wavelength accuracy          | ±2 nm   |
| Monochromator wavelength repeatability     | ±0.2 nm   |
| OD accuracy                                | <1% at 2.0 OD<br><3% at 3.0 OD  |
| OD linearity                               | <1% from 0 to 3.0 OD  |
| OD repeatability                           | <0.5% at 2.0 OD   |
| Stray light                                | 0.03% at 230 nm   |
| Reading speed (kinetic)                    | 96 well: 6 seconds<br>384 well: 11 seconds<br>1536 well: 25 seconds   |
| Fluorescence Intensity                     |   |
| Light source                               | Xenon flash   |
| Detector                                   | Dual top PMTs;<br>Single top PMT (option)<br>Low noise PMT (bottom filter system)<br>Red shifted PMT (top/bottom monochromator system |
| Wavelength selection                       | Quad monochromators (top/bottom)<br>Filters (top/bottom)  |
| Wavelength range                           | Monochromators: 250 – 850 nm<br>Filters (dual PMT): 200 – 850 nm  |
| Monochromator bandwidth                    | Variable; from 3 to 50 nm, in 1 nm increments   |
| Dynamic range                              | 7 decades   |

| Specifications are subject to change. Performance values represent the average |
|--|
| observed factory test values. See www.biotek.com for a complete list.          |

| Sensitivity (Fluorescein)  | Filters: 0.2 pM (4 amol/well, 384-well low vol plate) - top 1 pM (10 amol/well, 1536-well plate) - top 1 pM (0.1 fmol/well, 384-well plate) - bottom Quad Monochromator: 2 pM (40 amol/well, 384-well low vol plate) - top 2.5 pM (0.25 fmol/well, 384-well plate) - bottomm |
|----------------------------|--|
| Reading speed<br>(kinetic) | 96 well: 6 seconds<br>384 well: 11 seconds<br>1536 well: 25 seconds  |
| Luminescence               |  |
| Wavelength range           | 300 - 700 nm   |
| Dynamic range              | >6 decades   |
| Sensitivity                | 5 amol ATP (384-well low volume plate)   |
| Fluorescence Polariz       | ation  |
| Light source               | Xenon flash  |
| Detector                   | Dual PMT or single PMT (option)  |
| Wavelength selection       | Filters  |
| Wavelength range           | 280 - 850 nm   |
| Sensitivity                | 1 mP standard deviation at 1 nM fluorescein<br>(394-well low volume plate)<br>1.5 mP standard deviation at 1 nM fluorescein<br>(1536-well plate)   |
| Time-Resolved Fluor        | rescence   |
| Light source               | Xenon flash  |
| Detector                   | Dual PMT or single PMT (option)  |
| Wavelength selection       | Quad monochromators (top/bottom)<br>Filters (top/bottom)   |
| Wavelength range           | Monos: 250 – 850 nm<br>Filters (dual PMT): 200 – 850 nm  |
| Sensitivity                | Europium 40 fM (384-well low volume plate)<br>Europium 70 fM (1536-well plate)   |
| Alpha Detection            |  |
| Light source               | 100 mW 680 nm laser  |
| Detector                   | PMT  |
| Wavelength selection       | Filters (top)  |
| Sensitivity                | 100 amol bio-LCK-P (384-well low volume plate)   |
| Read speed                 | 96 well: 30 seconds<br>384 well: 1 minute 50 seconds<br>1536 well: 7 minutes 20 seconds  |
| Reagent Dispensers         |  |
| Number                     | 2 syringe pumps  |
| Supported labware          | 6- to 384-well plates, Petri dishes  |
| Dead volume                | 1.1 mL, 100 µL with back flush   |
| Dispense volume            | 5-1000 μL in 1 nm increments   |
| Physical Characteris       | tics   |
| Power                      | 250 Watts max.   |
| Dimensions                 | 15.4" W x 20.7" D x 16.1" H<br>(39 x 52.5 x 41 cm)   |
| Weight                     | 78 lbs (35Kg)  |
| Regulatory                 |  |
| Regulatory                 | CE and TUV marked. RoHS Compliant. Models for In Vitro Diagnostic use are available.   |

## Synergy™ H1 Multi-Mode Reader

Synergy™ H1 Multi-Mode Reader is equipped with both monochromator and filter optical systems. Synergy H1 provides flexibility and performance, at a very attractive price.

#### Flexibility at a Great Price

Synergy H1 is available in a monochromator-only configuration. Supporting top and bottom fluorescence, UV-visible absorbance and luminescence, it is the most cost-effective solution of its type on the market. Combined with the Take3™ Micro-Volume Plate for low volume 2 µL assays, it is the perfect instrument for lifescience research laboratories.

#### **Patented Hybrid Optical System**

Adding the optional filter module turns the Synergy H1 into an advanced Hybrid reader. This patented optical design is only available from BioTek. Monochromators provide ease-of-use and flexibility, while filters provide increased optical efficiency and sensitivity.

## Gas Controller for Live Cell Assays

An available Gas Controller for Synergy H1 allows control and monitoring of  $\mathrm{CO}_2$  and  $\mathrm{O}_2$  levels in the system. The Gas Controller, along with advanced temperature control to 45 °C and orbital shaking, create the ideal physiological environment needed

for assays using live cells. Live cell workflows can be automated by integrating Synergy H1 with BioSpa™ 8 Automated Incubator.

## Upgradable to Advanced Read Modes

When equipped with the optional filter module, Synergy H1 may be used for fluorescence polarization assays as well as Time-Resolved Fluorescence (TRF) and TR-FRET assays.

#### **Dual Reagent Injector Module**

For rapid, precise reagent injection in all plate types, Synergy H1 has an available dual reagent injector module, ideal for inject/read applications.

- Nucleic acid quantification
- Protein quantification
- Enzyme kinetics
- ▶ Biomarker quantification
- ELISAs
- Yeast kinetic analysis
- Genetic analysis
- Drug discovery
- Cell proliferation
- Cytotoxicity
- Drug absorption and metabolism
- ▶ Biologics drug discovery and development
- Food safety
- ▶ Environmental monitoring



Synergy H1 shown with optional Gas Controller module





## Specifications

| General                                    |  |
|--|--|
| Detection mode                             | UV-Vis absorbance<br>Fluorescence intensity<br>Luminescence<br>Fluorescence polarization<br>Time-resolved fluorescence |
| Read methods                               | Endpoint, kinetic, spectral scanning, well area scanning   |
| Microplate types                           | 1- to 384-well plates  |
| Other labware supported                    | Petri and cell culture dishes<br>Take3™ Micro-Volume Plates  |
| Temperature control                        | 4-Zone™ incubation to 45 °C with Condensation Control™   |
| Shaking                                    | Linear, orbital, double orbital  |
| Software                                   | Gen5 <sup>™</sup> Data Analysis Software<br>Gen5 Secure for 21 CFR Part 11 compliance (option)                         |
| Automation                                 | BioStack <sup>™</sup> and 3rd party automation compatible<br>BioSpa <sup>™</sup> 8 Automated Incubator compatible      |
| CO <sub>2</sub> and O <sub>2</sub> control | Range: 0 - 20% (CO <sub>2</sub> ); 1 - 19% (O <sub>2</sub> )   |
| Absorbance                                 |  |
| Light source                               | Xenon flash  |
| Detector                                   | Photodiode   |
| Wavelength selection                       | Monochromator  |
| Wavelength range                           | 230 to 999 nm, in 1 nm increments  |
| Monochromator<br>bandwidth                 | 4 nm (230-285 nm); 8 nm (>285 nm)  |
| Dynamic range                              | 0 - 4.0 OD   |
| Resolution                                 | 0.0001 OD  |
| Pathlength correction                      | Yes  |
| Monochromator wavelength accuracy          | ±2 nm  |
| Monochromator wavelength repeatability     | ±0.2 nm  |
| OD accuracy                                | <1% at 2.0 OD<br><3% at 3.0 OD   |
| OD linearity                               | <1% from 0 to 3.0 OD   |
| OD repeatability                           | <0.5% at 2.0 OD  |
| Stray light                                | 0.03% at 230 nm  |
| Reading speed (kinetic)                    | 96 wells: 11 seconds<br>384 wells: 22 seconds  |
| Fluorescence Intensity                     | 1  |
| Light source                               | Xenon flash  |
| Detector                                   | PMT for monochromator system<br>PMT for filter system  |
| Wavelength selection                       | Quad monochromators (top/bottom)<br>Filters (top)  |
| Wavelength range                           | Monochromators: 250 - 700 nm (850 nm option)<br>Filters: 200 - 700 nm (850 nm option)                                  |

| Monochromator<br>bandwidth   | Fixed, 16 nm  |
|------------------------------|---|
| Dynamic range                | 5 decades   |
| Sensitivity<br>(Fluorescein) | Filters:<br>0.25 pM (0.025 fmol/well, 384-well plate)   |
|                              | Quad Monochromator:<br>2.5 pM (0.25 fmol/well, 384-well plate) - top<br>4 pM (0.4 fmol/well, 384-well plate) - bottom |
| Reading speed (kinetic)      | 96 wells: 11 seconds<br>384 wells: 22 seconds   |
| Luminescence                 |   |
| Wavelength range             | 300 - 700 nm  |
| Dynamic range                | >6 decades  |
| Sensitivity                  | Monos: 20 amol ATP (flash)<br>Filters: 10 amol ATP (flash), 100 amol (glow)   |
| Fluorescence Polar           | ization   |
| Light source                 | Xenon flash   |
| Detector                     | PMT   |
| Wavelength selection         | Filters   |
| Wavelength range             | 280 - 700 nm (850 nm option)  |
| Sensitivity                  | 1.2 mP standard deviation at 1 nm fluorescein   |
| Time-Resolved Flu            | orescence   |
| Light source                 | Xenon flash   |
| Detector                     | PMT   |
| Wavelength selection         | Quad monochromators (secondary mode)<br>Filters (top)   |
| Wavelength range             | Filters: 200 - 700 nm (850 nm option)   |
| Sensitivity                  | Filters: Europium 40 fM (4 amol/well, 384-well plate)<br>Monos: Europium 1200 fM (120 amol/well, 384-well plate)      |
| Reagent Dispense             | rs  |
| Number                       | 2 syringe pumps   |
| Supported labware            | 6- to 384-well plates, Petri dishes   |
| Dead volume                  | 1.1 mL, 100 µL with back flush  |
| Dispense volume              | 5-1000 μL in 1 nm increments  |
| Dispense accuracy            | ±1 µL or 2%   |
| Dispense precision           | ≤2% at 50-200 µL  |
| Physical Character           | stics   |
| Power                        | 130 Watts max.  |
| Dimensions                   | 15.4"W 18.6"D 12.9"H<br>(39.1 x 47.2 x 32.8 cm)   |
| Weight                       | 50 lbs (22.5 kg)  |
| Regulatory                   |   |
| Regulatory                   | CE and TUV marked. RoHS Compliant. Models for In Vitro Diagnostic use are available.                                  |

### Synergy™ 2 Multi-Mode Reader

The Synergy™ 2 has been designed for life science research and drug discovery applications. It incorporates enhanced fluorescence, luminescence and absorbance optics for superior performance. Advanced read modes such as fluorescence polarization and time-resolved fluorescence are available as individual, upgradable modules and an available reagent injection system expands the applications range.

#### **Best Price/Performance Ratio**

The Synergy 2 incorporates dedicated, optimized optical paths for each detection mode using filters for fluorescence and a monochromator for absorbance. The result is excellent performance in all modes, at an attractive price.

## Sensitive Dichroic-based Fluorescence Optics

The Synergy 2 fluorescence optics are a step up from the Synergy HTX, incorporating dichroic mirrors, which decrease background noise, as well as a liquid-filled emission fiber that increases the system's light collection efficiency. The result is higher sensitivity for demanding assays.

#### Dedicated Luminescence Light Path

A dedicated liquid-filled light guide coupled with a low noise detector provides high-performance luminescence detection, on par with dedicated microplate luminometers.

Synergy 2 is DLReady™ certified by Promega to run their Dual-Luciferase® assay system.

#### Advanced, Modular Read Modes

In addition to the basic read modes available on the Synergy HTX, Synergy 2 offers fluorescence polarization, time-resolved fluorescence and Alpha detection modes, available as individual, upgradable modules.

- Nucleic acid quantification
- Protein quantification
- Enzyme kinetics
- ▶ Biomarker quantification
- ▶ ELISAs
- Genetic analysis
- Drug discovery
- ▶ Cell proliferation
- Cytotoxicity
- Drug absorption and metabolism
- ▶ Biologics drug discovery and development
- Food safety
- Biofuels research
- Environmental monitoring





### **Specifications**

| General                                |   |
|--|---|
| Detection mode                         | UV-Vis absorbance<br>Fluorescence intensity<br>Luminescence<br>Fluorescence polarization<br>Time-resolved fluorescence<br>Alpha             |
| Read methods                           | Endpoint, kinetic, spectral scanning, well area scanning  |
| Microplate types                       | 6- to 1536-well plates<br>1- to 384-well (luminescence)   |
| Other labware<br>supported             | Petri and cell culture dishes<br>Take3™ Micro-Volume Plates   |
| Temperature control                    | 4-Zone <sup>™</sup> incubation to 65 °C   |
| Shaking                                | Linear  |
| Software                               | Gen5 <sup>™</sup> Data Analysis Software<br>Gen5 Secure for 21 CFR Part 11 compliance (option)  |
| Automation                             | BioStack™ and 3rd party automation compatible   |
| Absorbance                             |   |
| Light source                           | Xenon flash   |
| Detector                               | Photodiode  |
| Wavelength selection                   | Monochromator   |
| Wavelength range                       | 200 to 999 nm, in 1 nm increments   |
| Monochromator<br>bandwidth             | 2.4 nm  |
| Dynamic range                          | 0 - 4.0 OD  |
| Resolution                             | 0.0001 OD   |
| Pathlength correction                  | Yes   |
| Monochromator wavelength accuracy      | ±2 nm   |
| Monochromator wavelength repeatability | ±0.2 nm   |
| OD accuracy                            | <1% at 2.0 OD<br><3% at 3.0 OD  |
| OD linearity                           | <1% from 0 to 3.0 OD  |
| OD repeatability                       | <0.5% at 2.0 OD   |
| Stray light                            | 0.03% at 230 nm   |
| Reading speed (kinetic)                | 96 wells: 11 seconds<br>384 wells: 22 seconds<br>1536 wells: 42 seconds   |
| Fluorescence Intensit                  | у   |
| Light source                           | Tungsten halogen<br>Xenon flash (option)  |
| Detector                               | PMT   |
| Wavelength selection                   | Filters/ dichroic mirrors   |
| Wavelength range                       | Tungsten lamp: 300 - 700 nm (850 nm option)<br>Xenon lamp: 200 - 700 nm (850 nm option)   |
| Dynamic range                          | >6 decades  |
| Sensitivity (Fluorescein)              | Top: 1 pM (0.2 fmol/well 96-well plate; 0.1 fmol/well 384-well plate) Bottom: 5 pM (1 fmol/well 96-well plate; 0.5fmol/well 384-well plate) |
| Reading speed (kinetic)                | 96 wells: 11 seconds<br>384 wells: 22 seconds<br>1536 wells: 42 seconds   |

| 300 - 700 nm  |
|---|
| >6 decades  |
| 10 amol ATP (flash)<br>100 amol ATP (glow)                              |
| ization   |
| Tungsten halogen<br>High energy DPR xenon flash (option)                |
| PMT   |
| Filters/dichroics   |
| 400 - 700 nm (320 - 850 nm option)                                      |
| 3 mP at 1 nM fluorescein  |
| orescence   |
| Xenon flash   |
| PMT   |
| Filters/dichroics   |
| Filters: 200 - 700 nm (850 nm option)                                   |
| Europium 60 fM (12 amol/well 96-well plate; 6 amol/well 384-well plate) |
|   |
| Tungsten halogen  |
| PMT   |
| Filters   |
| 100 amol of bio-LCK-P, 25 μL/ well in 384-well plate                    |
| 2 minute (96-well plate)  |
| rs  |
| All modes   |
| 2 syringe pumps   |
| 6- to 384-well plates   |
| 1.1 mL, 100 μL with back flush  |
| 5-1000 μL in 1 nm increments  |
| ±1 μL or 2%   |
| ≤2% at 50-200 μL  |
| istics  |
| 250 Watts max.  |
| 17"W x 17.5"D x 14.5"H  |
| (43.5 x 44.5 x 37.3 cm)   |
|   |
| (43.5 x 44.5 x 37.3 cm)   |
|   |

## Synergy™ HTX Multi-Mode Reader

The Synergy™ HTX is an entry-level, affordable and upgradeable multi-mode microplate reader. Available read modes include top and bottom fluorescence, UV-visible absorbance and luminescence detection.

Temperature control to 50 °C, shaking and advanced Gen5™ data analysis software are also included. A dual reagent injector module is available for all read modes and plate types.

#### Ideal for Basic Research Applications

The Synergy HTX is the ideal instrument for nucleic acid and protein quantification, enzyme assays, biomarker quantification and ELISA assays, as well as cell-based assays (gene expression, cellular growth, cytotoxicity).

#### AlphaScreen®/ AlphaLISA®

AlphaScreen and AlphaLISA assays can be performed on Synergy HTX with excellent results. Alpha-capable configurations add assay versatility to basic research requirements.

## Sensitive Filter-based Fluorescence

Two excitation and two emission filters are included with the reader, and can be used for top and bottom reading. Bottom reading is usually recommended when working with adherent cells, as it often provides better signal-to-background ratios. Top reading is typically best for assays where the fluorescence signal comes from the solution.

#### Flexible Monochromator-based Absorbance

All Synergy readers use monochromators for absorbance detection. This provides unlimited wavelength selection from the low UV to the near infrared, in 1 nm steps, and enables spectral scanning.

### Low-noise Luminescence Detection

The Synergy HTX can automate glow and flash luminescence assays, thanks to its optional dual reagent injector module. Typical assays include ATP quantification as well as luciferase gene expression assays.

- ► AlphaScreen®/AlphaLISA®
- Nucleic acid quantification
- Protein quantification
- Enzyme kinetics
- Biomarker quantification
- ▶ ELISAs
- ► Genetic analysis
- ▶ Cell proliferation
- Cytotoxicity
- Drug absorption and metabolism
- Food safety
- Environmental monitoring





## **Specifications**

| =                                      |   |
|--|---|
| General                                |   |
| Detection mode                         | UV-Vis absorbance<br>Fluorescence intensity<br>Luminescence<br>Time-resolved fluorescence (secondary mode)<br>Alpha |
| Read methods                           | Endpoint, kinetic, spectral scanning, well area scanning  |
| Microplate types                       | 6- to 384-well plates   |
| Other labware<br>supported             | PCR plates, Petri and cell culture dishes<br>Take3™ Micro-Volume Plates   |
| Temperature control                    | 4-Zone™ incubation to 50 °C with Condensation Control™  |
| Shaking                                | Linear, orbital   |
| Software                               | Gen5™ Data Analysis Software<br>Gen5 Secure for 21 CFR Part 11 compliance (option)                                  |
| Automation                             | BioStack™ and 3rd party automation compatible   |
| Absorbance                             |   |
| Light source                           | Xenon flash   |
| Detector                               | Photodiode  |
| Wavelength selection                   | Monochromator   |
| Wavelength range                       | 200 to 999 nm, in 1 nm increments   |
| Monochromator<br>bandwidth             | 2.4 nm  |
| Dynamic range                          | 0 - 4.0 OD  |
| Resolution                             | 0.0001 OD   |
| Pathlength correction                  | Yes   |
| Monochromator wavelength accuracy      | ±2 nm   |
| Monochromator wavelength repeatability | ±0.2 nm   |
| OD accuracy                            | <1% at 2.0 OD<br><3% at 3.0 OD  |
| OD linearity                           | <1% from 0 to 3.0 OD  |
| OD repeatability                       | <0.5% at 2.0 OD   |
| Stray light                            | 0.03% at 230 nm   |
| Reading speed (kinetic)                | 96 wells: 14 seconds<br>384 wells: 26 seconds   |
| Fluorescence Intensit                  | у   |
| Light source                           | Tungsten halogen<br>Xenon flash (option)  |
| Detector                               | PMT   |
| Wavelength selection                   | Filters   |
| Wavelength range                       | 300 - 700 nm (200 - 850 nm option)  |
| Dynamic range                          | >6 decades  |
| Sensitivity (Fluorescein)              | Top and Bottom: 5 pM (1 fmol/well 96-well plate)  |
| Reading speed (kinetic)                | 96 wells: 31 seconds<br>384 wells: 80 seconds   |
|  |   |

| Wavelength range          | 300 - 700 nm  |
|---------------------------|---|
| Dynamic range             | >6 decades  |
| Sensitivity               | 10 amol ATP (flash) - Lum and Abs / Lum configurations<br>30 amol ATP (flash) - Multi-mode configurations |
| Time-Resolved Flu         | prescence   |
| Light source              | Xenon flash   |
| Detector                  | PMT   |
| Wavelength selection      | Monochromator   |
| Alpha Detection           |   |
| Light source              | Tungsten halogen  |
| Detector                  | PMT   |
| Wavelength selection      | Filters   |
| Sensitivity               | 300 amol of bio-LCK-P, 25 μL/ well in 384-well plate  |
| Read speed                | 2 minute (96-well plate)  |
| Reagent Dispenser         | rs .  |
| Supported detection modes | All modes   |
| Number                    | 2 syringe pumps   |
| Supported labware         | 6- to 384-well plates   |
| Dead volume               | 1.1 mL, 100 µL with back flush  |
| Dispense volume           | 5-1000 μL in 1 nm increments  |
| Dispense accuracy         | ±1 µL or 2%   |
| Dispense precision        | ≤2% at 50-200 μL  |
| Physical Characteri       | stics   |
| Power                     | 100-240 VAC, 50/60 Hz   |
| Dimensions                | 6"W x 15"D x 10"H<br>(40.6 x 38 x 25.4 cm)  |
| Weight                    | 40 lbs (18 kg)  |
| Regulatory                |   |
| Regulatory                |   |

## Multi-Mode Reader Comparison Chart

Which Multi-Mode Reader is right for you?

|   | Synergy™ Neo2   | Cytation™  |
|---|---|--|
| Key Features  |   |  |
| Monochromator-based<br>UV-visible absorbance        | •   | •  |
| Fluorescence top/bottom                             | •   | •  |
| Luminescence  | •   | •  |
| Filtered luminescence                               | •   | •  |
| Injectors   | •   | •  |
| TRF & TR-FRET                                       | •   | •  |
| Fluorescence polarization                           | •   | •  |
| Standard AlphaLISA®/AlphaScreen®                    |   |  |
| Laser AlphaLISA/AlphaScreen                         | •   | Cytation 5 only  |
| Hybrid Technology™                                  | •   | •  |
| Dual PMT read head                                  | •   |  |
| Performance Specifications                          |   |  |
| Fluorescein typical – top                           | 2 pM (monos) / 0.2 pM (filters)                       | 2.5 pM (monos) / 0.25 pM (filters)   |
| Fluorescein typical – bottom                        | 1 pM (filters)  | 4 pM (monos)   |
| ATP typical – flash luminescence                    | 5 amol  | 10 amol  |
| Polarization typical – 1 nM Fluorescein             | 1 mP ST   | 1.2 mP SD  |
| Europium typical                                    | 40 fM   | 40 fM (filters)  |
| AlphaScreen typical - LCK peptide                   | 100 amol  | 100 amol (Cytation 5)  |
| Fastest read speed<br>96-/384-well plates (seconds) | 6/11  | 11/22  |
| General Specifications                              |   |  |
| Microplate types                                    | 6 to 1536   | Mono: 6 to 384<br>Filter and imaging: 6-1536   |
| Gas Controller compatible                           | •   | •  |
| BioSpa™ 8 Automated Incubator compatible            | •   | •  |
| Automation ready/BioStack™ compatible               | •   | •  |
| Dual reagent injector compatible                    | •   | •  |
| Barcode reader option                               | •   |  |
| Take3™ Micro-Volume Plate compatible                | •   | •  |
| Temperature control system                          | to 65 °C  | To 45° C (Cytation 3)<br>To 65° C (Cytation 5)   |
| Condensation Control™                               | •   | •  |
| Filter capacity                                     | Up to 6 filter sets                                   | 2 filter sets  |
| Fluorescence bandwidth                              | Filter dependent<br>Mono: variable from 3 nm to 50 nm | Filter dependent<br>Mono: variable from<br>9 nm to 50 nm<br>(Cytation 5)<br>Fixed 16 nm (Cytation 3) |

| Synergy H1                         | Synergy 2        | Synergy HTX                                     |  |
|------------------------------------|------------------|---|--|
|                                    |                  |   |  |
| •                                  | •                | •   |  |
| •                                  | •                | •   |  |
| •                                  | •                | •   |  |
| •                                  | •                | •   |  |
| •                                  | •                | •   |  |
| •                                  | •                | (secondary mode)                                |  |
| •                                  | •                |   |  |
|                                    | •                | •   |  |
|                                    |                  |   |  |
| •                                  |                  |   |  |
|                                    |                  |   |  |
| 2.5 pM (monos) / 0.25 pM (filters) | 1 pM             | 5 pM  |  |
| 4 pM (monos)                       | 5 pM             | 5 pM  |  |
| 10 amol                            | 10 amol          | 30 amol/10 amol (luminescence on configuration) |  |
| 1.2 mP SD                          | 3 mP             | n/a   |  |
| 40 fM (filters)                    | 60 fM            | n/a   |  |
| n/a                                | 100 amol         | 300 amol  |  |
| 11/22                              | 11/22            | 14/26   |  |
|                                    |                  |   |  |
| 6 to 384                           | 6 to 1536        | 6 to 384  |  |
| •                                  |                  |   |  |
| •                                  |                  |   |  |
| •                                  | •                | •   |  |
| •                                  | •                | •   |  |
|                                    |                  |   |  |
| • 45 9C                            | • /F 9C          | •<br>• FO 9C                                    |  |
| to 45 °C                           | to 65 ℃          | to 50 °C  |  |
| •                                  |                  | •   |  |
| 2 filter sets                      | 4 filter sets    | 4 filter sets                                   |  |
| Filter dependent<br>Mono 16 nm     | Filter dependent | Filter dependent                                |  |

### FLx800™ Fluorescence Reader

The compact FLx800™ fluorescence reader provides high performance in 6- to 384-well microplates at an attractive price. Options include top and bottom detection, temperature control and reagent injector.

#### **Great Price/Performance Ratio**

The FLx800 uses top and bottom bifurcated quartz fibers to ensure strong sample excitation and efficient collection of the emitted signal. This reader combines sensitivity, convenience and ease of use all at a great price.

#### **Bottom Reading and FRET Detection**

The bottom detection system uses a large 5 mm diameter quartz fiber optimized for cell-based assays in 96-well plates and smaller densities. The reader may be equipped with up to 4 filter sets and may be used to run cell-based FRET assays.

#### Sensitive Luminescence Detection

The FLx800's detector provides very high sensitivity when running luminescent assays. ATP or luciferase can be quantified down to very low concentrations using the reader's photon integration mode.

#### Fluorescent Ion Channel Assays

A syringe pump injector is available as an option to automate fluorescent ion-channel assays. This system is used to inject a trigger reagent that induces a fast change in fluorescent signal. The FLx800 kinetically monitors the signal just after injection.

### Typical Applications:

- Nucleic acid quantification
- Protein quantification
- Enzyme kinetics
- ▶ Genetic analysis by fluorescence
- ► Cellular analysis by fluorescence



### **Specifications**

| General                   |  |  |  |
|---------------------------|--|--|--|
| Detection mode            | Fluorescence intensity<br>Luminescence   |  |  |
| Read method               | Endpoint, kinetic, well area scanning (under computer control)                                 |  |  |
| Microplate types          | 6- to 384-well plates  |  |  |
| Temperature control       | 4 °C above ambient to 50 °C (I models)   |  |  |
| Shaking                   | Linear (I models)  |  |  |
| Software                  | Gen5 <sup>™</sup> Data Analysis Software<br>Gen5 Secure for 21 CFR Part 11 compliance (option) |  |  |
| Fluorescence Intensi      | ity  |  |  |
| Light source              | Tungsten halogen   |  |  |
| Detector                  | PMT  |  |  |
| Wavelength selection      | Filters  |  |  |
| Wavelength range          | 300 - 700 nm (850 nm option)   |  |  |
| Dynamic range             | 5 decades  |  |  |
| Sensitivity               | Fluorescein 5 pM (1 fmol/well 96-well plate)   |  |  |
| Luminescence              |  |  |  |
| Wavelength range          | 300 - 700 nm   |  |  |
| Dynamic range             | 5 decades  |  |  |
| Sensitivity               | 100 amol ATP (flash)   |  |  |
| Reagent Dispensers        |  |  |  |
| Supported detection modes | Fluorescence   |  |  |
| Number                    | 1 syringe pump   |  |  |
| Dead volume               | 1.1 mL, 100 μL   |  |  |
| Dispense volume           | 5-1000 μL in 1 nm increments   |  |  |
| Physical Characteris      | tics   |  |  |
| Power                     | 100 - 240 Volts AC. 50/60 Hz.  |  |  |
| Dimensions                | 15"W x 16"D x 9"H<br>(38.1 x 40.64 x 22.89 cm)   |  |  |
| Weight                    | 30 lbs (13.6 kg  |  |  |
| Regulatory                |  |  |  |
| Regulatory                | CE and TUV marked. RoHS compliant. Models for In Vitro Diagnostic use are available.           |  |  |

## PowerWave™ HT Microplate Spectrophotometer

The PowerWave™ HT is a high throughput, robot friendly microplate spectrophotometer with a very small footprint, ideal for integration into automated systems. Wavelength selection in 1 nm increments, temperature control and superior performance up to 4.0 OD add to its appeal for a variety of assay needs.

#### High Speed, Higher Throughput

In automation and high throughput, timing is everything... With 8 reading channels, the PowerWave HT can read a 96-well plate in 5 seconds.

## Low Stray Light Monochromator Optics

PowerWave HT's monochromator optics pre-select the measurement wavelength before light goes through the sample. This results in very low stray light reaching the detector...with the added benefit of excellent performance even at high optical densities.

#### BioStack<sup>™</sup> Compatible for Benchtop Automation

When walkaway benchtop automation is required, the PowerWave HT, coupled with BioStack, provides a compact system for rapid processing of up to 50 plates at a time.

#### Gen5™ Control = Assay Flexibility

Gen5 Data Analysis Software not only allows easy control of all the functionality of the PowerWave HT, it also supports a vast number of applications in absorbance. Quick export to Microsoft® Excel® or use Gen5's powerful data analysis tools to make quick work of the most complex assays.

#### Typical Applications:

- ► Enzyme kinetics
- ▶ ELISAs
- Genetic analysis by colorimetry
- Cellular analysis by colorimetry
- Cell proliferation



### **Specifications**

| pecifications                                |   |  |  |
|--|---|--|--|
| General                                      |   |  |  |
| Detection mode                               | Absorbance  |  |  |
| Read method                                  | Endpoint, kinetic   |  |  |
| Microplate types                             | 96- and 384-well plates   |  |  |
| Temperature control                          | 4-Zone™ incubation to 50 °C   |  |  |
| Shaking                                      | Linear  |  |  |
| Software                                     | Gen5™ Data Analysis Software<br>Gen5 Secure for 21 CFR Part 11 compliance (option)  |  |  |
| Automation                                   | BioStack <sup>™</sup> and 3rd party automation compatible                           |  |  |
| Absorbance                                   |   |  |  |
| Light source                                 | Xenon flash   |  |  |
| Detector                                     | Photodiodes   |  |  |
| Wavelength selection                         | Monochromator   |  |  |
| Wavelength range                             | 200 to 999 nm, in 1 nm increments   |  |  |
| Monochromator<br>bandwidth                   | 5 nm  |  |  |
| Dynamic range                                | 0 - 4.0 OD  |  |  |
| Resolution                                   | 0.001 OD  |  |  |
| Pathlength correction                        | Yes   |  |  |
| Monochromator wavelength accuracy            | ±2 nm   |  |  |
| Monochromator<br>wavelength<br>repeatability | ±0.2 nm   |  |  |
| OD accuracy                                  | 1% + 0.01 OD  |  |  |
| OD linearity                                 | ±1%   |  |  |
| OD repeatability                             | 0.5% + 0.005 OD   |  |  |
| Stray light                                  | 0.03% at 230 nm   |  |  |
| Reading speed<br>(kinetic)                   | 96 wells: 5 seconds<br>384 wells: 11 seconds  |  |  |
| Physical Characteris                         | stics   |  |  |
| Power  | 100-240 Volts AC. 50/60 Hz  |  |  |
| Dimensions                                   | 8.5"W x 16"D x 8.5"H<br>(21.6 x 40.6 x 21.6 cm)                                     |  |  |
| Weight                                       |   |  |  |
| Regulatory                                   |   |  |  |
| Regulatory                                   | CE and TUV marked. ROHS compliant.<br>In Vitro Diagnostic use models are available. |  |  |

### Epoch™ 2 Microplate Spectrophotometer

Epoch™ 2 is a compact monochromator-based microplate spectrophotometer for 6- to 384-well microplates, cuvettes and 2 µL measurements. Epoch 2 features a 10" color touchscreen interface with easy to navigate controls, and full onboard Gen 5™ software for data collection, powerful analysis and flexible export and report options. Incubation, shaking and robot compatibility are standard features.

#### **UV-Vis Measurements**

Epoch 2's monochromator-based optics offers wavelength selection from 200 nm to 999 nm – for applications from nucleic acid quantification to ELISA, without using filters. Epoch 2 can measure up to 48 2  $\mu$ L samples in the unique Take3 Micro-Volume plates for rapid direct quantification. An optional cuvette port provides

quick 1 cm measurements, making Epoch 2 a versatile spectrophotometer for multiple applications.

#### Touch. Run. Done.

Designed for easy-to-use, yet powerful functionality, Epoch 2 features a 10" color touchscreen interface, WiFi, Bluetooth and USB connectivity and flash drive storage. It's a self-contained computer, in a space and cost saving design, configurable for the laboratory's needs today and in the future.

#### Full Gen5 Data Analysis Software

With Epoch 2, "onboard software" doesn't mean "limited software". Complete reader control, protocol design, data analysis and export/report functions are at your fingertips. For applications in microplates, cuvettes or

Take3 plates, Gen5 offers the same intuitive navigation and full capability as an external computer. With Gen5 on the Epoch 2 – there's no need for a dedicated computer – it's all built-in!

#### Advanced 4-Zone™ Incubation

Epoch 2 features BioTek's 4-Zone natural convection incubator up to 65 °C with minimal variation across the plate – ideal for a wide range of temperature-sensitive assays. Epoch 2's unique Condensation Control™, solves the common problem of condensation build up on plate lids during incubated kinetic runs. Epoch 2 can be integrated with BioSpa™ 8 Automated Incubator for unattended automation.

- ELISA
- Enzyme kinetics
- ▶ Nucleic acid and protein quantification
- Cell proliferation
- Cytotoxicity
- Spectral scanning
- Reactive oxygen species
- ► Food safety and quality
- ▶ Bacterial identification
- ► Total protein determination
- Nucleic acid purity assessment





### **Specifications**

| General                                |   |  |
|--|---|--|
| Detection mode                         | Absorbance  |  |
| Read method                            | Endpoint, kinetic, well area scanning   |  |
| Microplate types                       | 6- to 384-well plates   |  |
| Other labware supported                | Take3™ Micro-Volume plates  |  |
| Temperature control                    | 4-Zone™ incubation to 65 °C   |  |
| Shaking                                | Linear, orbital, double-orbital   |  |
| Software                               | Gen5™ Data Analysis Software<br>Gen5 Secure for 21 CFR Part 11 compliance (option)  |  |
| Automation                             | BioStack™ and 3rd party automation compatible<br>BioSpa™ 8 Automated Incubator compatible   |  |
| Absorbance                             |   |  |
| Light source                           | Xenon flash   |  |
| Detector                               | Photodiode  |  |
| Wavelength selection                   | Monochromator   |  |
| Wavelength range                       | 200 to 999 nm, in 1 nm increments   |  |
| Monochromator bandwidth                | 2.9 nm  |  |
| Dynamic range                          | 0 - 4.0 OD  |  |
| Resolution                             | 0.0001 OD   |  |
| Pathlength correction                  | Yes   |  |
| Monochromator wavelength accuracy      | ±2 nm   |  |
| Monochromator wavelength repeatability | ±0.2 nm   |  |
| OD accuracy                            | 0 to 2 OD: +1% +0.010 OD<br>2 to 2.5 OD: +3% +0.010 OD  |  |
| OD linearity                           | 0 to 2.0 OD ± 1% ± 0.010<br>2.0 to 2.5 OD ± 3% ± 0.010  |  |
| OD repeatability                       | 0 to 2.0 OD ± 1% ± 0.005<br>2.0 to 2.5 OD ± 3% ± 0.005  |  |
| Stray light                            | 0.03% at 230 nm   |  |
| Reading speed (kinetic)                | 96 wells: 8 seconds<br>384 wells: 14 seconds  |  |
| Physical Characteristics               |   |  |
| Power                                  | 100-240 Volts AC. 50/60 Hz  |  |
| Dimensions                             | With touchscreen: 15.5" D x 12.5" W x 13" H (39.3 x 32 x 33 cm) Without touchscreen: 15.5" D x 12.5" W x 8" H (39.3 x 32 x 20.3 cm) |  |
| Weight                                 | With touchscreen: 25 lbs (11.3 kg) Without touchscreen: 20 lbs (9.1 kg)   |  |
| Regulatory                             |   |  |
| Regulatory                             | CE and TUV marked, RoHS compliant. Models for In Vitro Diagnostic use are available.  |  |

### Epoch™ Microplate Spectrophotometer

Epoch™ is a monochromator-based microplate spectro-photometer that offers superior functionality for the life science laboratory at an accessible price. Controlled by the powerful, yet easy-to-use Gen5™ Data Analysis Software, Epoch is designed to be the new lab workhorse for a wide variety of applications. For walk-away automation, an optional BioStack™ compatible Epoch is available.

#### 200 nm to 999 nm Wavelength Range

The monochromator-based optical system in Epoch allows any wavelength selection between 200

and 999 nm in 1 nm increments. No filters required! From low UV nucleic acid measurements to standard ELISA assays, Epoch is ideally suited to the life science laboratory where application flexibility is required.

#### 6- to 384-well Microplate Reading

Epoch's optical and mechanical systems are designed to provide optimal measurements in a variety of microplates. The area scanning capability provides multiple measurements across larger diameter wells, resulting in more meaningful data analysis.

## Take3™ Micro-Volume Plate Compatible

When sample size matters, as in critical nucleic acid and protein quantification, the Take3 plate provides up to sixteen 2  $\mu$ L measurements – without needing to dilute important samples.

## **Endpoint, Kinetic, Spectral Scanning**

There's no need to buy expensive instrumentation to perform a variety of absorbance measurements. Epoch, driven by Gen5 Data Analysis Software, is the ultimate high-value system with maximum assay flexibility.

- Nucleic acid quantification
- ▶ Protein quantification
- ▶ 260/280 and 260/230 purity measurements
- ► ELISA
- Enzyme kinetics
- Cytotoxicity
- Cell proliferation
- ▶ Micro-volume assays with Take3 plate





### **Specifications**

| General                                |  |  |
|--|--|--|
| Detection mode                         | Absorbance   |  |
| Read method                            | Endpoint, kinetic, well area scanning  |  |
| Microplate types                       | 6- to 384-well plates  |  |
| Other labware supported                | Take3™ Micro-Volume plates   |  |
| Software                               | Gen5™ Data Analysis Software<br>Gen5 Secure for 21 CFR Part 11 compliance (option) |  |
| Automation                             | BioStack™ and 3rd party automation compatible ("R"model)                           |  |
| Absorbance                             |  |  |
| Light source                           | Xenon flash  |  |
| Detector                               | Photodiode   |  |
| Wavelength selection                   | Monochromator  |  |
| Wavelength range                       | 200 to 999 nm, in 1 nm increments  |  |
| Monochromator bandwidth                | 5 nm   |  |
| Dynamic range                          | 0 - 4.0 OD   |  |
| Resolution                             | 0.0001 OD  |  |
| Pathlength correction                  | yes  |  |
| Monochromator wavelength accuracy      | ±2 nm  |  |
| Monochromator wavelength repeatability | ±0.2 nm  |  |
| OD accuracy                            | 0 to 2 OD: +1% +0.010 OD<br>2 to 2.5 OD: +3% +0.010 OD                             |  |
| OD linearity                           | 0 to 2.0 OD ± 1% ± 0.010<br>2.0 to 2.5 OD ± 3% ± 0.010                             |  |
| OD repeatability                       | 0 to 2.0 OD ± 1% ± 0.005<br>2.0 to 2.5 OD ± 3% ± 0.005                             |  |
| Reading speed (kinetic)                | 96 wells: 15 seconds<br>384 wells: 31 seconds                                      |  |
| Physical Characteristics               |  |  |
| Power                                  | 100-240 Volts AC. 50/60 Hz   |  |
| Dimensions                             | 12" W x 12.5" D x 7.7" H<br>(30.5 cm x 31.8 cm x 19.6 cm)                          |  |
| Weight                                 | <15 lbs (6.8 kg)   |  |
| Regulatory                             |  |  |
| Regulatory                             | CE and TUV marked. ROHS compliant. In Vitro Diagnostic use models are available.   |  |

### ELx808™ Absorbance Reader

Excellent optical performance and superior incubation are among the top features of this multi-channel reader. The ELx808™ is suitable for a wide array of applications, from endpoint ELISAs to kinetic cell growth studies.

#### **4-Zone™ Temperature Control**

For temperature sensitive assays, there is no better incubation system in this microplate reader class than the ELx808. The natural convection heating is software controlled for consistency and performance over time.

#### **Fast Measurement**

The ELx808 can collect kinetic data in intervals as short as 6 seconds, for the most demanding assays. Gen5™ Data Analysis Software provides multiple kinetic and end point data analysis options for a variety of applications.

#### **Superior Optical Performance**

The ELx808 can accommodate up to six absorbance filters, and its optical channels are staggered to prevent crosstalk between wells. The reference channel eliminates channel-to-channel variation. This unique design gives the ELx808 its proven optical performance.

#### Typical Applications:

- ► ELISA
- Enzyme kinetics
- ▶ Endotoxin assays
- Cell growth studies
- Cytotoxicity
- Protein assays



### **Specifications**

| General                 |  |  |  |
|-------------------------|--|--|--|
| Detection mode          | Absorbance   |  |  |
| Read method             | Endpoint, kinetic, linear scanning                                       |  |  |
| Microplate types        | 96-well plates   |  |  |
| Temperature control     | 4-Zone™ incubation to 50 °C<br>±0.2 °C at 37 °C                          |  |  |
| Software                | Gen5™ Reader Control Software<br>Gen5™ Data Analysis Software (optional) |  |  |
| Absorbance              |  |  |  |
| Light source            | Tungsten halogen   |  |  |
| Detector                | Photodiode   |  |  |
| Wavelength selection    | Filters  |  |  |
| Wavelength range        | 380 to 900 nm<br>340 to 900 nm (ELx808IU)                                |  |  |
| Dynamic range           | 0 - 4.0 OD   |  |  |
| Resolution              | 0.001 OD   |  |  |
| Pathlength correction   | No   |  |  |
| OD accuracy             | <1% at 2.5 OD<br><2% at 3.5 OD   |  |  |
| OD linearity            | <1% at 2.5 OD  |  |  |
| OD repeatability        | <0.5% at 2.5 OD<br><1.5% at 3.5 OD                                       |  |  |
| Reading speed (kinetic) | 96 wells: 8 seconds  |  |  |
| Physical Characterist   | tics   |  |  |
| Power                   | 100-240 Volts AC. 50/60 Hz   |  |  |
| Dimensions              | 15.5"W x 16"D x 8.75"H<br>(39.4 x 40.6 x 22.2 cm)                        |  |  |
| Weight                  | 30 lbs (13.6 kg)   |  |  |
| Regulatory              |  |  |  |
| Regulatory              | For In Vitro Diagnostic use. CE and TUV marked, RoHS Compliant.          |  |  |



### ELx800™ Absorbance Reader

The ELx800™ is a compact, robust microplate reader ideally suited for applications within the clinical and life science research laboratories.

#### Reliable And Robust Design

These characteristics are the reason there have been more than 25,000 ELx800s placed in laboratories around the globe. There simply isn't a more reliable reader with proven performance than the ELx800.

#### Gen5™ Software Expands Versatility

Under computer control by Gen5 Data Analysis Software, ELx800 applications are expanded to include kinetic and well area scanning measurements. Data analysis and reporting/exporting features in Gen5 are completely customizable to suit your laboratory's requirements.

#### High Performance, High Value

The ELx800 comes with four filters (user-selectable wavelengths), with a five filter capacity. In addition to standard 96-well microplate reading, the ELx800 offers 6-, 12-, 24-, 48- and optional 384-well microplate reading, to fit a variety of assay needs.

### **Typical Applications:**

- ► ELISA
- Protein assays
- Cytotoxicity



### **Specifications**

| General                 |  |  |
|-------------------------|--|--|
| Detection mode          | Absorbance   |  |
| Read method             | Endpoint. Kinetic, well area scanning (under computer control)                                   |  |
| Microplate types        | 6- to 384-well plates  |  |
| Software                | Gen5 <sup>™</sup> Reader Control Software<br>Gen5 <sup>™</sup> Data Analysis Software (optional) |  |
| Absorbance              |  |  |
| Light source            | Tungsten halogen   |  |
| Detector                | Photodiode   |  |
| Wavelength selection    | Filters  |  |
| Wavelength range        | 400 to 750 nm<br>340 to 750 nm (UV option)   |  |
| Dynamic range           | 0 - 3.0 OD   |  |
| Resolution              | 0.001 OD   |  |
| Pathlength correction   | No   |  |
| OD accuracy             | <1% at 2.0 OD  |  |
| OD linearity            | <1% at 2.0 OD<br><3% at 3.0 OD   |  |
| OD repeatability        | <0.5% at 2.0 OD  |  |
| Reading speed (kinetic) | 96 wells: 30 seconds   |  |
| Physical Characterist   | tics   |  |
| Power                   | 100-240 Volts AC. 50/60 Hz   |  |
| Dimensions              | 15"W x 16.5"D x 7"H<br>(38.1 x 41.9 x 17.8 cm)   |  |
| Weight                  | 18.5 lbs (8 kg)  |  |
| Regulatory              |  |  |
| Regulatory              | For In Vitro Diagnostic use. CE and TUV marked, RoHS Compliant.                                  |  |



## Absorbance Reader Comparison Chart

### Which Absorbance Reader is right for you?

|   | PowerWave™ HT          | Epoch™2  | Epoch™   |  |
|---|------------------------|--|--|--|
| Key Features                                |                        |  |  |  |
| Wavelength selection                        | Monochromator<br>based | Monochromator<br>based                                     | Monochromator<br>based                                     |  |
| Wavelength range (nm)                       | 200 - 999              | 200 - 999  | 200 - 999  |  |
| Microplate types                            | 96 and 384             | 6 to 384   | 6 to 384   |  |
| Absorbance range                            | 0 - 4.0                | 0 - 4.0  | 0 - 4.0  |  |
| Temperature control                         | to 50 °C               | to 65 °C   |  |  |
| Shaking                                     | Linear                 | Linear, orbital, double-orbital                            |  |  |
| Cuvette measurement                         | Cuvette adapter        | Cuvette port (optional),<br>Take3 or cuvette adapter       | Take3 or cuvette adapter                                   |  |
| Filter capacity                             | n/a                    | n/a  | n/a  |  |
| Automation ready/BioStack™<br>compatible    | •                      | •  | ("R" configuration)  |  |
| BioSpa™ 8 Automated Incubator compatible    |                        | •  |  |  |
| Gen5 <sup>™</sup> Software version included | Gen5                   | Gen5   | Gen5   |  |
| Take3™ Micro-Volume Plate compatible        |                        | •  | •  |  |
| Fastest read speed: 96 wells (seconds)      | 5                      | 8  | 15   |  |
| Typical Performance                         |                        |  |  |  |
| OD accuracy                                 | 1% +0.01 OD            | 0 to 2.0 OD: ±1% ±0.010 OD<br>2.0 to 2.5 OD: ±3% ±0.010 OD | 0 to 2.0 OD: +1% +0.010 OD<br>2.0 to 2.5 OD: +3% +0.010 OD |  |
| OD linearity                                | ±1%                    | 0 to 2.0 OD: ±1% ±0.010 OD<br>2.0 to 2.5 OD: ±3% ±0.010 OD | 0 to 2.0 OD: +1% +0.010 OD<br>2.0 to 2.5 OD: +3% +0.010 OD |  |
| OD repeatability                            | 0.5% ±0.005 OD         | 0 to 2.0 OD: ±1% ±0.005 OD<br>2.0 to 2.5 OD: ±3% ±0.005 OD | 0 to 2.0 OD: +1% +0.005 OD<br>2.0 to 2.5 OD: +3% +0.005 OD |  |
| Resolution                                  | 0.001 OD               | 0.0001 OD  | 0.0001 OD  |  |

| ELx808 IU                          | ELx808™                            | ELx800™                        | ELx800 UV                      | ELx800 NB                   |
|------------------------------------|------------------------------------|--------------------------------|--------------------------------|-----------------------------|
|                                    |                                    |                                |                                |                             |
| Filter-based                       | Filter-based                       | Filter-based                   | Filter-based                   | Filter-based                |
|                                    |                                    |                                |                                |                             |
| 340 - 900                          | 380 - 900                          | 400 - 750                      | 340 - 750                      | 400 - 750                   |
| 96                                 | 96                                 | 6 to 96                        | 6 to 96                        | 6 to 384                    |
| 0 - 4.0                            | 0 - 4.0                            | 0 - 3.0                        | 0 - 3.0                        | 0 - 3.0                     |
| to 50 °C                           |                                    |                                |                                |                             |
| Linear                             | Linear                             |                                |                                |                             |
|                                    |                                    |                                |                                |                             |
|                                    |                                    |                                |                                |                             |
| 6                                  | 6                                  | 5                              | 5                              | 5                           |
|                                    |                                    |                                |                                |                             |
|                                    |                                    |                                |                                |                             |
|                                    |                                    |                                |                                |                             |
|                                    |                                    |                                |                                |                             |
| Gen5RC                             | Gen5RC                             | Gen5RC                         | Gen5RC                         | Gen5RC                      |
|                                    |                                    |                                |                                |                             |
| 8                                  | 8                                  | 30                             | 30                             | 30                          |
| 0                                  | 0                                  | 30                             | 30                             | 30                          |
|                                    |                                    |                                |                                |                             |
| <1% at 2.5 OD                      | <1% at 2.5 OD                      | <1% at 2.0 OD                  | <1% at 2.0 OD                  | <1% at 2.0 OD               |
| <2% at 3.5 OD                      | <2% at 3.5 OD                      |                                |                                |                             |
| <1% at 2.5 OD                      | <1% at 2.5 OD                      | <1% at 2.0 OD<br><3% at 3.0 OD | <1% at 2.0 OD<br><3% at 3.0 OD | <1% at 2.0 OD <3% at 3.0 OD |
| 0.50/05.05                         | .0.50′ . 0.5.05                    |                                |                                |                             |
| <0.5% at 2.5 OD<br><1.5% at 3.5 OD | <0.5% at 2.5 OD<br><1.5% at 3.5 OD | <0.5% at 2.0 OD                | <0.5% at 2.0 OD                | <0.5% at 2.0 OD             |
|                                    |                                    | 0.001.00                       | 0.001.00                       | 0.001.00                    |
| 0.001 OD                           | 0.001 OD                           | 0.001 OD                       | 0.001 OD                       | 0.001 OD                    |

### Gen5™ Data Analysis Software

Gen5™ Data Analysis Software incorporates over 30 years of experience and user feedback into outstanding microplate reader software. Gen5 is a unique combination of power and ease-of-use that drives productivity and saves time. Use Gen5 to control BioTek's readers and export data, or as a fully integrated processing tool.

#### **Beginner-friendly Software**

Gen5 is built around logical laboratory workflows to read microplates and produce/analyze data. In Gen5, you simply click "Read Now" and follow the prompts. At the end of the read, answer the question: "Do you want to export to Excel?" With Gen5 you don't have to spend hours figuring out how to get things done.

#### **Powerful Functionality**

Gen5 comes with powerful built-in tools such as 4-P and 5-P curve fits with or without weighting, parallel-line analysis, advanced kinetic analysis, and much more. The software has been specifically designed to analyze matrices of data that are difficult to process in Microsoft® Excel® spreadsheets. Special attention has been placed on result presentation so complex data can be displayed in a clean, colorful way to facilitate data interpretation.

#### Up-to-date Web-based Sample Files


A searchable library of sample files with data is available on the BioTek web site. Existing files are kept up to date and new files are added on a regular basis. Gen5 users can upload files to share their experiments and sample data with other Gen5 users.

#### Gen5 Secure: You Are In Control


Gen5 Secure is for you if secured data storage, user group management or 21 CFR Part 11 compliance are a core part of your requirements. It includes extra features such as 25 licenses per copy, quality control trending module with Levey-Jennings charts and automatic email notification on trigger events.

#### All-in-One Solution


Control any BioTek reader, automate the entire process and produce publication-ready reports with one integrated, powerful software. It doesn't get more efficient than this.



Parallel line analysis and  $EC_{50}$  determinations



High-resolution 99 x 99 area scan



Multiple window views of plate data and results

## $\textbf{Gen5}^{\scriptscriptstyle\mathsf{TM}}\ \textbf{Comparison}\ \textbf{Chart}$

|   | Gen5 RC   | Gen5      | Gen5 Image+ | Gen5 Secure | Gen5 Secure Image+ |
|---|-----------|-----------|-------------|-------------|--------------------|
| Readers Supported   |           |           |             |             |                    |
| All Cytation and Synergy readers; Epoch 2 (external),<br>Epoch, PowerWave HT, FLx800, ELx808, ELx800  | •         | •         | •           | •           | •                  |
| Read Modes and Methods  |           |           |             |             |                    |
| Absorbance, fluorescence, luminescence, time-resolved fluorescence, fluorescence polarization, Alpha  | •         | •         | •           | •           | •                  |
| End point, kinetic, spectral scanning, linear scanning, Take3 interface   | •         | •         | •           | •           | •                  |
| Single and Multi-Mode Data Analysis   |           |           |             |             |                    |
| Qualitative, quantitative, kinetic, spectral analysis, custom transformations, EC50, parallel line analysis, $Z^\prime$ calculation, validation and cutoff formulae |           | •         | •           | •           | •                  |
| Imaging Modes and Methods   |           |           |             |             |                    |
| Fluorescence, brightfield, phase contrast, color brightfield  | •         | •         | •           | •           | •                  |
| Single color, multi-color, montage, Z-stacking  | •         | •         | •           | •           | •                  |
| Image Processing and Analysis   |           |           |             |             |                    |
| Cell counting   |           | •         | •           | •           | •                  |
| Subpopulation analysis  |           |           | •           |             | •                  |
| Hit-picking   |           |           | •           |             | •                  |
| Z-projection  |           |           | •           |             | •                  |
| Digital phase contrast  |           |           | •           |             | •                  |
| Image stitching   |           |           | •           |             | •                  |
| Security / 21 CFR Part 11   |           |           |             |             |                    |
| User groups, single sign-on (SSO) option  |           |           |             | •           | •                  |
| Secure database data storage, audit trails  |           |           |             | •           | •                  |
| Electronic signature, email notification  |           |           |             | •           | •                  |
| Validation Tools  |           |           |             |             |                    |
| Gen5 Validation Package   | available | available | available   | available   | available          |

### Take3™ Micro-Volume Plate

Quickly quantify ultra-low volume samples of DNA, RNA and protein. Measure up to 48 samples with volumes as low as 2 µL without dilution. Take3™ can be used to measure a standard cuvette or patented BioCells™ for quick 1 cm measurements. Low volume, higher throughput is available with the Take3 Trio.

## Compatible with Most BioTek Detection Systems

Epoch™, Synergy™ and Cytation™ reader functionality can easily reach into the micro-volume range using the Take3 plate. Measure multiple 2 μL samples, cuvettes or BioCells. Adding the Take3 plate to a BioTek detection system creates an incredibly versatile workstation for a variety of applications.

## Unique Robust Construction and Easy Maintenance

The anodized aluminum base construction, precision crafted slides and hydrophobic sample surfaces make

#### **Typical Applications:**

- Micro-volume DNA, RNA and protein quantification
- Micro-volume fluorescence measurements in Synergy and Cytation readers
- ► Fluorescent dye incorporation measurements
- Spectral scanning in micro-volume, cuvette or BioCell



pipetting simple and cleanup effortless. For routine cleaning of the sample surfaces, a laboratory wipe is all that's needed. If a slide becomes damaged, replacement is easy – no need to return the Take3 to the factory for repair or calibration.

### Gen5 Take3 Module: Automated DNA, RNA and Protein Quantification

It couldn't be easier to get multiple (up to 48) nucleic acid or protein sample results. Gen5's Take3 module includes pre-programmed protocols with immediate results output including spectral scans and purity ratios. There's no need for complicated configuration or calculation.

### **Specifications**

|                      | Take3         | Take3 Trio    |
|----------------------|---------------|---------------|
| 2 μL sample capacity | 16            | 48            |
| Detection limit      | 2 ng/μL dsDNA | 2 ng/μL dsDNA |
| BioCell capacity     | 2             | 2             |
| Cuvette capacity     | 1             | n/a           |

Specifications are subject to change.

### Reader Accessories

BioTek offers a wide range of accessories to help increase productivity, expand your plate reader's capabilities, and maintain the performance of your BioTek microplate reader system. See our web site for a complete listing of available accessories.



#### **Dual Reagent Injector Module**

Automate inject/read assays such as flash luminescence assays (ATP, luciferase) and fluorescent ion channel assays on all Synergy and Cytation readers.



#### **Gas Controller**

The Gas Controller module for the Synergy H1, Synergy Neo2 and Cytation allows full control over CO<sub>2</sub> and O<sub>2</sub> concentrations to modulate the environment for microplate-based live cell assays.



#### **Gen5 Secure Software**

Upgrade to Gen5 Secure for 21 CFR Part 11 compliance, user management features, data encryption and much more.



**Filters and Mirrors** 

A full range of standard and custom filters and dichroic mirrors are available for applications from the low UV to the near infrared.



#### Instrument Qualification

See the Compliance Section on pages 62-63 for details about BioTek's product qualification tools and services.



**BioStack™ Microplate Stacker** 

Automate routine processes with this compact stacker. BioStack is also compatible with BioTek's liquid handling instruments.



 $\textbf{BioSpa}^{\text{\tiny{M}}} \textbf{ 8 } \textbf{ Automated Incubator}$ 


BioSpa 8 is an automated incubator linking BioTek readers or imagers together with washers and dispensers for full workflow automation of up to 8 microplates.





"I have used the MultiFlo FX Microplate Dispenser since August 2014 for my routine research experiments and assays such as ELISAs and other biochemical kinetic/endpoint reactions. The dispenser is indispensable in modern laboratories... It is unbelievably compact, saving the much needed workbench space. The software and colour touch-screen user interface have been well designed for the average user - with fast and simple protocol programming and operation making it the simplest machine to operate...The Dispenser has been very effective in my research, particularly on the quality of my experiments. "

(MultiFlo FX Multi-Mode Dispenser)



BioTek offers a range of compact and affordable solutions for your laboratory's specific liquid handling needs. BioTek is world renowned for manufacturing the most reliable and versatile microplate washers on the market. From basic ELISA to sensitive cell washing to bead washing (including Luminex® xMAP® technology), the EL406™, 405™ Touch, ELx50™ and MultiFlo™ FX are configured with many options to meet a myriad of assay requirements. For unattended automation of live cell and other assays, the 405, EL406 and MultiFlo FX integrate with the BioSpa™ 8 Automated Incubator. From milliliters down to one microliter, our reagent dispensers offer simple, repeatable and precise liquid delivery throughout their defined volume range. Single, 8and 12-channel transfer tools are available, along with bulk reagent dispensers, to meet varied liquid handling requirements.

### EL406™ Washer Dispenser

The EL406™ Combination
Washer Dispenser is the only
instrument on the market
offering fast microplate washing
together with BioTek's unique
Parallel Dispense™ technologies
for optimized liquid handling
processes.

# Unattended Automation of ELISAs and Cell-based Assays

The EL406 integrates 96-, 384- and 1536-well microplate washing with three dispensers in one compact instrument. Now you can simply press a button and walk away, or automate an entire batch by adding a BioStack™ Microplate Stacker. For unattended automation the EL406 can be integrated to BioSpa™ 8 Automated Incubator. Add a BioTek imager or reader and entire workflows can be automated.

# Patented Dual-Action™ Manifold and Ultrasonic Advantage™

The EL406 incorporates BioTek's Dual-Action manifold for thorough yet gentle washing of loosely adherent cell layers, and Ultrasonic Advantage for automated wash manifold maintenance.

### **Parallel Dispense Technologies**

The EL406 eliminates the need to choose a dispensing technology by offering both peristaltic and syringe pumps on a single platform.

### Fast and Efficient Biomagnetic Separation and Vacuum Filtration

The EL406 automates full microplate washing of magnetic microspheres used in an increasing number of multiplex assays. Developed in conjunction with Luminex® xMAP® technology leaders, BioTek's separation modules incorporate high energy neodymium iron boron magnets for speed and efficiency. An available vacuum filtration module makes the EL406 also well suited for polystyrene beads and filtration-to-waste processes.

- ▶ ELISA automation
- ▶ MSD assay automation
- ▶ High content screening immunocytochemistry
- ▶ Cell-based assays
- ▶ FLIPR® Ca<sup>2+</sup> flux
- Magnetic bead assay automation
- ▶ Polystyrene bead assay automation
- Drug transport assays
- Automated cell washing, fixing and staining for cellular imaging
- SiLA compliant integration (with LHC software)





| General                     |  |  |
|-----------------------------|--|--|
| Microplate types            | 96-, 384-, 1536-well<br>Low profile and standard height<br>Solid and filter bottom (option)  |  |
| Onboard software            | Create, edit or run multiple protocols   |  |
| Software<br>(pc control)    | LHC Software (option) LHC Secure for 21 CFR Part 11 compliance (option) SiLA Compliant driver (option)   |  |
| Separation                  | Biomagnetic separation, vacuum filtration (option)   |  |
| Shaking                     | Programmable up to 60 minutes<br>Slow, medium, fast or variable  |  |
| Soaking                     | Programmable up to 60 minutes  |  |
| Ultrasonic<br>Advantage™    | Yes (standard on most configurations)  |  |
| Automation                  | BioStack <sup>™</sup> and 3rd party automation compatible<br>BioSpa <sup>™</sup> 8 Automated Incubator compatible  |  |
| Washing                     |  |  |
| Manifold types              | 96-well washing: 96-tube manifold<br>96- and 384-well washing: 96-tube Dual-Action manifold<br>384-well washing (fast): 192-tube Dual-Action manifold<br>1536-well washing: Two 32-tube dispense manifolds, 316<br>SS tubes or sapphire jeweled 316 SS tubes |  |
| Volume range                | 3- 3,000 μL/well, in 1 μL increments   |  |
| Wash cycles                 | 1-250  |  |
| Buffer/reagent selection    | Auto switching module for up to 4 buffers (option)   |  |
| Supply bottle               | 4 L or 10 L (optional)   |  |
| Dispense accuracy           | ±3%  |  |
| Dispense precision          | <3% CV (model dependent)   |  |
| Residual volume             | <2 μL/well   |  |
| Wash speed                  | 96 wells, 300 μL/well, 96-tube manifold: 13 seconds 384 wells, 100 μL/well, 192-tube manifold: 17 seconds 1536 wells, 10 μL/wells, two 32-tube manifolds: 36 seconds   |  |
| Flow rates                  | High flow to low flow<br>Optimized rates for cell assays   |  |
| Sterilization               | Chemical   |  |
| Vacuum range for filtration | 0 to -380 mmHg   |  |
| Dispensing - peri           | staltic pump (multi-channel)   |  |
| Manifold types              | 8-tip (1 x 8) cassette with plastic, 316 stainless steel or sapphire jeweled 316 stainless steel tips  |  |
| Dispense speed              | 96 wells, 10 $\mu$ L /well: 8 seconds 384 wells, 5 $\mu$ L /well: 12 seconds 1536 wells, 1 $\mu$ L /well: 27 seconds   |  |
| Volume range                | 500 nL - 3,000 μL/well, selectable in 1 μL increments  |  |
| Flow rates                  | User programmable rates from high to low<br>Optimized rates for cell assays  |  |

| Dispense<br>performance                           | 1 μL: Recommended Volume Range: 1 - 50 μL Dispense Accuracy: ±5% at 1 μL Dispense Precision: ≤5% CV at 1 μL ≤10% CV at 500 nL Minimum Prime Volume: 1.20 mL 5 μL cassette: recommended range: 5 - 2,500 μL Dispense Accuracy: ±2.0% at 5 μL Dispense Precision: ≤2.5% CV at 5 μL Minimum Prime Volume: 4.23 mL 10 μL cassette: recommended range: 10 - 3,000 μL Dispense Accuracy: ±2.0% at 10 μL Dispense Precision: 2.0 CV at 10 μL Minimum Prime Volume: 7.36 mL |  |
|---|---|--|
| Recommended<br>cassette replace-<br>ment interval | 1 μL Cassette: 1,000 384-well microplates at 5 μL well 5 μL Cassette: 1,000 96-well microplates at 50 μL well 10 μL Cassette: 1,000 96-well microplates at 100 μL well  |  |
| Sterilization                                     | Autoclave, chemical   |  |
| Dispensing - syringe pump (multi-channel)         |   |  |
| Manifold types                                    | 96-well dispensing: One 16-tube (2 x 8) manifold - 316 stainless steel tubes 96-/384-well dispensing: Two 16-tube (1 x 16) manifolds - 316 stainless steel tubes 1536-well dispensing: Two 32-tube (1 x 32) manifolds - sapphire jeweled 316 stainless steel or 316 stainless steel tubes   |  |
| Dispensing speed                                  | 20 μL /well, 96 wells, 1 x 16 tubes: 5 seconds<br>20 μL /well, 384 wells, 1 x 16 tubes: 14 seconds<br>3 μL /well 1536 wells, 2 x 32 tubes: 7 seconds  |  |
| Volume range                                      | 3 - 3,000 μL/well, selectable in 1 μL increments<br>Minimum Prime Volume: 12 mL   |  |
| Flow rates  | User programmable rates from high to low  |  |
| Dispense<br>accuracy                              | ±1 μL at 5 μL<br>±1 μL at 20 μL<br>±1% at 100 μL  |  |
| Dispense<br>precision                             | ≤5% CV at 5 μL<br>≤2.5% CV at 20 μL<br>≤1% CV at 100 μL   |  |
| Supply bottle                                     | 1 L or 2 L  |  |
| Sterilization                                     | Chemical, autoclavable optoin   |  |
| Physical Charac                                   | teristics   |  |
| Power   | 100 - 240 Volts AC. 50/60 Hz  |  |
| Dimensions  | 16.5" W x 18" D x 12.5" H (42 x 46 x 32 cm)   |  |
| Weight  | 32 lbs (14.5 kg)  |  |
| Regulatory  |   |  |
| Regulatory  | CE and TUV marked. ROHS compliant. In Vitro Diagnostic use models are available.  |  |
|   |   |  |

### 405™ Touch Washer

BioTek's 405™ Touch Microplate Washer takes plate washing to the next level with an enhanced user interface, increased convenience, assay applications and automated maintenance features.

#### Industry Leading, User Pleasing

The 405 Touch Microplate Washer incorporates all the features and functionality of the prior ELx405 models, and improves accessibility through its touchscreen and extensive onboard software.

96- and 384-well microplate based wash procedures are only 'two touches' away with the easy-to-use interface. Additionally, two USB flash drives provide convenient file transfer, storage and operation.

A context sensitive Help System and several instructional videos are also included.

#### The Standard for Automation

The 405 Microplate Washer makes quick work of any washing assay, and is especially well suited for integration into automated systems, where the wash process is controlled remotely. The 405 can be integrated with the BioSpa™ 8 Automated Incubator for unattended automation of many common processes.

### **Cell and Bead Assays**

The 405 is available in various models for optimized performance with the most sensitive and rigorous assay requirements. When the protocol calls for washing loosely adherent cells, the Select model is fine-tuned with angled dispense tubes, extra low flow rates and unique X, Y and Z positioning. Magnetic and polystyrene bead washing are effectively accomplished with the 405.

### Verify™ Technology and Automated Ultrasonic Cleaning

BioTek's new patented Verify technology runs an automated QC check for manifold tube blockage, and visually reports any failures. Patented Ultrasonic Advantage™ can then be used to automatically clean the manifolds. Together, these features make the 405 a self-checking, self-maintaining microplate washer!

## Applications in Deep Well Washing

The ELx405 Select Deep Well washes 96- and 384-well plates up to 50 mm tall, and is also compatible with standard height plates without any hardware or software changes. This versatile system is optimal for labs working in deep well blocks and standard plates.

- ELISAs
- MSD assays
- ► HCS immune cytochemistry
- ► FLIPR® Ca<sup>2+</sup> flux
- ► Cell-based assays
- Magnetic and polystyrene bead assays
- Gene expression assays
- Cytokine assays
- ▶ ELISPOT assays
- Plasmid DNA purification
- Serum/plasma sample preparation
- Cell signaling phospho flow setup for flow cytometry
- SiLA compliant integration (with LHC software)







| General                          | General   |  |  |
|----------------------------------|---|--|--|
| Microplate types                 | 96- and 384-well Low profile and standard height Solid and filter bottom (option) - Filter pore sizes from 0.45 µm to 1.2 µm                              |  |  |
| Onboard software                 | Create, edit or run multiple protocols  |  |  |
| Software (pc control)            | LHC Software LHC Secure for 21 CFR Part 11 compliance (option) SiLA Compliant driver (option)   |  |  |
| Separation                       | Biomagnetic separation, vacuum filtration (optional)  |  |  |
| Shaking                          | Programmable up to 60 minutes<br>Slow, medium, fast or variable   |  |  |
| Soaking                          | Programmable up to 60 minutes   |  |  |
| Ultrasonic Advantage™ available  | Yes   |  |  |
| Verify™ clog detection available | Yes   |  |  |
| Automation                       | BioStack <sup>™</sup> and 3rd party automation compatible<br>BioSpa <sup>™</sup> 8 Automated Incubator compatible   |  |  |
| Washing                          |   |  |  |
| Manifold types                   | 96-tube manifold for 96-well washing<br>96-tube Dual-Actionmanifold for 96- & 384-well washing<br>192-tube Dual-Action manifold for fast 384-well washing |  |  |
| Volume range                     | 25-3,000 μL/well, in 1 μL increments  |  |  |
| Wash cycles                      | 1-250   |  |  |
| Buffer/reagent selection         | Auto switching (internal) for up to 4 buffers (option)  |  |  |
| Supply bottle                    | 4 L or 10 L (optional)  |  |  |
| Dispense precision               | <3% CV: 300 μL/well (96-well washing)<br><4% CV: 80 μL/well (384-well washing)  |  |  |
| Residual volume                  | < 2 μL/well (96- & 384-well plates)<br>96-tube manifold for 96 wells; 192-tube for 384 wells  |  |  |
| Wash speed                       | 96-wells, 300 µL/well, 3 cycles; ≤30 seconds 384-wells,100 µL/well, 3 cycles: ≤80 seconds 384-wells, 400 µL/well, 1 cycle: ≤20 seconds                    |  |  |
| Flow rates                       | High flow to low flow<br>Optimized rates for cell assays  |  |  |
| Sterilization                    | Chemical  |  |  |
| Vacuum range for filtration      | -38 mmHg to -506 mmHg   |  |  |
| Ultrasonic Advantage™ available  | Option  |  |  |
| Verify™ clog detection available | Option  |  |  |
| Physical Characteristics         |   |  |  |
| Power                            | 100 - 240 Volts AC. 50/60 Hz  |  |  |
| Dimensions                       | 14" W x 17" D x 10" H (35.6 x 43.2 x 25.4 cm)   |  |  |
| Weight                           | With internal buffer switching - 36 lbs (16.5 kg)   |  |  |
| Regulatory                       |   |  |  |
| Regulatory                       | CE and TUV marked. ROHS compliant. In Vitro Diagnostic use models are available.  |  |  |
|                                  |   |  |  |

### ELx50™ Strip Washer

The ELx50™ Microplate Strip Washer's compact footprint conceals a powerhouse of washing capabilities unsurpassed in its class. Exceptional dispense precision and evacuation efficiency can be utilized for both 96- and 384-well strip and plate washing.

### Multiple Washers in One

The ELx50 is a 3-in-1 solution for 96-well plate formats automating the wash steps of ELISAs, magnetic bead assays and polystyrene bead assays. The ELx50 can also be equipped with BioTek's Dual-Action™ manifold allowing independent control of dispense and aspiration manifolds

in both 96- and 384-well formats. As a welcome upgrade from manual processing, the ELx50 provides an all-inclusive wash solution offering consistent performance and unattended operation.

### Syringe Drive Fluid Delivery

As a self-contained and programmable washer, the ELx50 allows for complete control of precise fluidic delivery from the gentle dripping of a simple squeeze bottle to the full force of pressure delivery systems. Comprehensive onboard software makes creating protocols quick and intuitive.

### **Automated Liquid Level Sensing**

Liquid Level Alert™ allows the convenience of continuous monitoring for both supply and waste bottles. At the beginning and end of a wash protocol, the liquid level is verified to ensure an adequate buffer remains to complete a wash. Sufficient storage capacity in the waste bottle is also verified.

- ELISAs
- ► Cell-based assays
- Magnetic bead assays
- Polystyrene bead assays
- ELISPOT assays
- Multiplex assays







| General                     |   |  |
|-----------------------------|---|--|
| Microplate types            | 96- and 384-well  |  |
| Onboard software            | Create, edit or run multiple protocols  |  |
| Separation                  | Biomagnetic separation, vacuum filtration (model dependent)   |  |
| Shaking                     | Programmable up to 60 minutes<br>Selectable intensity (15 - 19 cycles/second)   |  |
| Soaking                     | Programmable up to 60 minutes   |  |
| Washing                     |   |  |
| Manifold types              | 8-tube manifold for 96-well washing<br>12-tube manifold for 96-well washing<br>16-tube Dual-Action manifold for 96- & 384-well washing  |  |
| Volume range                | 25-3,000 μL/well, in 1 μL increments  |  |
| Wash cycles                 | 1-10  |  |
| Buffer/reagent selection    | Auto switching module for up to 3 buffers (option)  |  |
| Supply bottle               | 2 L   |  |
| Dispense precision          | ≤3% CV  |  |
| Residual volume             | Solid bottom plates: $\leq 2~\mu$ L/well Filter bottom plates: Average increase weight of plate $\leq 1.2~g$ after dispensing 300 $\mu$ L to 0.45 $\mu$ m plate, filtration 30 seconds, low vacuum, blotted |  |
| Wash speed                  | 96-wells, 8-tube manifold, >300 μL/well: <130 seconds   |  |
| Flow rates                  | High flow to low flow   |  |
| Sterilization               | Chemical  |  |
| Vacuum range for filtration | -91 mmHg to -313 mmHg   |  |
| Physical Characteristics    |   |  |
| Power                       | 100 - 240 Volts AC. 50/60 Hz  |  |
| Dimensions                  | 14"W x 16"D x 6.5"H (35.6 x 40.6 x 16.5 cm)   |  |
| Weight                      | 22 lbs (9.8 kg)   |  |
| Regulatory                  |   |  |
| Regulatory                  | CE and TUV marked. ROHS compliant. In Vitro Diagnostic use models are available.  |  |

# Washer Comparison Chart

| Which Washer is right for you?  | EL406™  | MultiFlo™ FX                 |
|---|---|------------------------------|
| Key Features  |   |                              |
| ELISA   | •   | •                            |
| Cell-based assays   | •   | •                            |
| Magnetic bead assays  | •   | •                            |
| Polystyrene bead assays   | •   |                              |
| Filtration-to-waste processes   | •   |                              |
| Touchscreen user interface  |   | •                            |
| USB ports for protocol transfer   |   | •                            |
| Performance Specifications  |   |                              |
| Washing speed: 3 aspirate/dispense cycles,<br>96-well solid bottom plate, 300 µL/well | ≤30 sec   | ≤130 sec                     |
| Dispense precision  | ≤3% CV  | ≤3% CV                       |
| Residual volume: solid bottom plate   | ≤2 µL/well  | ≤2 µL/well                   |
| Vacuum filtration: 1.2 µm 96-well plate   | 0 mmHg to -380 mmHg   |                              |
| General Specifications  |   |                              |
| Microplate types  | 96, 384 and 1536  | 6, 12, 24, 48, 96 and 384    |
| Low profile and standard height   | •   | •                            |
| Solid and filter bottom   | •   |                              |
| Deep well   |   |                              |
| Manifold  |   |                              |
| 6-, 12-, 24-, 48-well washing   |   | Custom manifolds available   |
| 96-well washing   | 96-tube (8x12)  |                              |
| 96-/384-well washing  | Dual-Action 96-tube (8x12)                                  | Dual-Action 8-tube (1x8)     |
| 384-well washing  | Dual-Action 192-tube (16x12)                                |                              |
| 1536-well washing   | Dispense: Two 32-tube (1x32)<br>Aspiration: 128-tube (4x32) |                              |
| Ultrasonic Advantage™   | •   |                              |
| Verify <sup>™</sup> technology  |   |                              |
| BioSpa™ 8 Automated Incubator compatible  | •   | •                            |
| Automation ready/BioStack™ compatible   | •   | •                            |
| Automatic buffer switching  | Up to 4   |                              |
| Flow rates  | High flow rates to low cell wash rates                      | High flow rates to low rates |
| Volume range  | 3 - 3000 μL/well  | 25 - 30,000 μL/well          |
| Microplate shaking  | •   | •                            |
| Fluid and waste detection   | •   | (optional)                   |
| Flow and vacuum detection   | •   | * 1                          |
| Overflow protection   | •   | •                            |
| Pre-programmed maintenance routines   | •   | •                            |
| Onboard software included   | •   | •                            |
| Liquid Handling Control™ Software compatible  | •   | •                            |



| • • • • • • • • • • • • • • • • • • •  |
|--|
| • • • • • • • • • • • • • • • • • • •  |
| • • • • • • • • • • • • • • • • • • •  |
|  |
|  |
| ≤30 sec       ≤30 sec       ≤130 sec         ≤3% CV       ≤3% CV       ≤3% CV         ≤2 μL/well       ≤2 μL/well       ≤2 μL/well   |
| ≤30 sec       ≤30 sec       ≤130 sec         ≤3% CV       ≤3% CV       ≤3% CV         ≤2 μL/well       ≤2 μL/well       ≤2 μL/well   |
| ≤3% CV ≤3% CV ≤3% CV ≤3% CV ≤3% CV ≤2 μL/well ≤2 μL/well ≤2 μL/well  |
| ≤3% CV ≤3% CV ≤3% CV ≤3% CV ≤3% CV ≤2 μL/well ≤2 μL/well ≤2 μL/well  |
| ≤2 μL/well ≤2 μL/well ≤2 μL/well ≤2 μL/well  |
|  |
| -38 mmHg to -506 mmHg  |
|  |
|  |
| 96 and 384 96 and 384 96 and 384 96 and 384  |
| • • • •  |
| •  |
| •  |
|  |
|  |
| 96-tube (8x12) 96-tube (8x12) 8-tube (1x8) or 12-tube (1x12)   |
| Dual-Action 96-tube (8x12)  Dual-Action 96-tube (8x12)  Dual-Action 96-tube (8x12)  Dual-Action 16-tube (1x16)   |
| Dual-Action 192-tube (16x12) Dual-Action 192-tube (16x12)  |
|  |
|  |
| •  |
| • • •  |
| •  |
| •  |
| Up to 4  |
| High flow rates to lowHigh flow rates to lowHigh flow rates to lowHigh flow rates to lowHigh flow rates to lowcell wash ratescell wash ratescell wash ratescell wash rates |
| 25 - 3000 μL/well 25 - 3000 μL/well 50 - 3000 μL/well 25 - 3000 μL/well  |
| •  |
| • • • •  |
| •  |
| • • • •  |
| •  |
| •  |
| •  |

## MultiFlo™ FX Multi-Mode Dispenser

The new MultiFlo™ FX Multi-Mode Dispenser offers modular configurations of up to four independent dispensers and a microplate washer in one compact platform. A color touchscreen interface makes programming quick and easy.

### Multi-Mode Dispensing

MultiFlo FX becomes a versatile multi-mode dispenser with the addition of either the RAD™ technology for random access dispensing to 6- to 384-well plates or a wash module for 6- to 384-well plate washing. Fast, intuitive programming and operation are via the color touchscreen user interface. A MultiFlo FX configured with either RAD™ technology or the wash module replaces up to five liquid handlers.

### Parallel Dispense™ Technologies

Offering BioTek's unique combination of peristaltic and

microprocessor controlled syringe pump dispensing, the MultiFlo FX enables users to choose which is best for a specific reagent. While peristaltic pumps offer low prime volumes and backflush capabilities, BioTek's syringe drives are programand-forget solutions that never require recalibration. Automated dispensing with walk-away confidence.

### Modular and Upgradable

The MultiFlo FX is configurable and upgradable from dispense or wash only, to a combined dispense and wash combination, with 1-to-4 reagent dispensing with peristaltic and/or syringe driven precision dispensers, or the addition of a RAD module for single channel and high volume dispensing. Purchase the modules required now, and upgrade in the future as assay needs change. With its compact footprint and base height of less than 8 inches, the MultiFlo FX

comfortably fits on any lab bench or robotic system. The MultiFlo FX is easily integrated with a BioStack™ Microplate Stacker for walk-away automation. For complete live cell workflow automation, MultiFlo FX can be integrated with the BioSpa™ 8 Automated Incubator.

# Versatile Applications - Liquid Handling

A wide array of plate type settings accommodates 6- to 1536-well plate formats up to 50 mm high for dispensing. Volumes from 500 nL µL to 3 mL are dispensed with accuracy and precision. The wash module works with 6- to 384-well plates in standard, half and deep well, and with cluster or mini-tubes. Automate cell-based assays by integrating MultiFlo FX with the new BioStack™ 4 to handle lidded plates with speed and ease

- ► Cell-based assays
- Primary and secondary screening assays
- Dispense/wash protocol automation
- Compound storage
- ▶ Genomics and proteomics research
- Magnetic bead assays
- ▶ ELISAs
- Multiplex assays
- Automated cell washing, fixing and staining for cellular imaging



| Missaulata turasa                             | Dianamaina (nasi numa and numinas).  |  |
|---|--|--|
| Microplate types                              | Dispensing (peri pump and syringe): 96-, 384- and 1536-well standard, deep and PCR plates 6 12-, 24- and 48-well plates (dispense tip configurable) Washing: 96 384-well standard plates; 6-, 12-, 24-, and 48-well plates (with compatible manifold)  |  |
| Onboard software                              | Create, edit or run multiple protocols   |  |
| Software<br>(pc control)                      | LHC Software<br>LHC Secure for 21 CFR Part 11 compliance (option)<br>SiLA Compliant driver (option)  |  |
| Shaking                                       | Programmable up to 60 minutes<br>Slow, medium, fast or variable  |  |
| Soaking                                       | Programmable up to 60 minutes  |  |
| Automation                                    | BioStack™ and 3rd party automation compatible<br>BioSpa™ 8 Automated Incubator compatible  |  |
| Washing                                       |  |  |
| Manifold types                                | 96- and 384-well washing: 8-tube manifold<br>Custom manifolds available for 6-, 12-, 24-, 48-well washing  |  |
| Volume range                                  | 20-30,000 μL/well  |  |
| Wash cycles                                   | 1-10   |  |
| Supply bottle                                 | 2 L  |  |
| Dispense accuracy                             | ±3%  |  |
| Dispense precision                            | $\leq$ 3% CV (96-/384-well plates; 300 µL/well) $\leq$ 5% CV (6-well plates; 5560 µL/well)   |  |
| Residual volume                               | < 2 μL/well, 300 μL dispense, 0.1% Tween   |  |
| Wash speed                                    | 96-wells, 8-tube manifold, >300 µL/well: <130 seconds  |  |
| Flow rates                                    | 140 to 422 µL/well/second  |  |
| Sterilization                                 | Chemical   |  |
| Dispensing - peristaltic pump (multi-channel) |  |  |
| Manifold types                                | $1\times8$ - sapphire jeweled 316 SS, 316 SS or polypropylene tips, with 1, 5 or 10 $\mu L$ tubing   |  |
| Fluid delivery                                | 1 or 2 peristaltic pumps   |  |
| Dispense speed                                | 96 wells, 5 µL cass, 10 µL/ well: 3 seconds<br>96 wells, 5 µL cass, 20 µL/ well: 3.5 seconds<br>384 wells, 5 µL cass, 5 µL/ well: 6.5 seconds<br>384 wells, 10 µL cass, 10 µL/ well: 8 seconds<br>384 wells, 1 µL cass, 1 µL/ well: 6 seconds<br>1536 wells, 1 µL cass, 1 µL/ well: 21 seconds   |  |
| Volume range                                  | 500 nL - 3,000 $\mu$ L/well, selectable in 1 $\mu$ L increments  |  |
| Flow rates                                    | User programmable rates from high to low   |  |
| Dispense<br>performance                       | 1 μL: Recommended Volume Range: 1 - 50 μL Dispense Accuracy: ±5% at 1 μL Dispense Precision: ≤5% CV at 1 μL ≤10% CV at 500 nL Minimum Prime Volume: 0.78 - 1.20 mL 5 μL: Recommended Volume Range: 5 - 2,500 μL Dispense Accuracy: ±2.0% at 5 μL Dispense Precision: ≤2.5% CV at 5 μL Minimum Prime Volume: 2.75 - 4.23 mL 10 μL: Recommended Volume Range: 10 - 3,000 μL Dispense Accuracy: ± 2.0% at 10 μL Dispense Precision: ≤-2.0% CV at 10 μL Minimum Prime Volume: 4.79 - 7.36 mL |  |
| Recommended                                   | 1 μL Cassette: 1,000 384-well microplates at 5 μL well<br>5 μL Cassette: 1,000 96-well microplates at 50 μL well   |  |
| cassette replacement<br>interval              | 5 μL Cassette: 1,000 96-well microplates at 50 μL well<br>10 μL Cassette: 1,000 96-well microplates at 100 μL well   |  |

| Dispensing - syring           | ge pump (multi-channel)  |  |
|-------------------------------|--|--|
| Manifold types                | 96- and 384-well dispensing:<br>One 16-tube (2 $\times$ 8) manifold - 316 SS tubes<br>Two 16-tube (1 $\times$ 16) manifolds - 316 SS tubes<br>1536-well dispensing:<br>Two 32-tube (1 $\times$ 32) manifolds - sapphire jeweled 316 SS or 316 SS tubes<br>6- to 48-well dispensing: custom autoclavable manifolds available  |  |
| Fluid delivery                | Two positive displacement syringe drives   |  |
| Dispensing speed              | 20 $\mu$ L /well, 96 wells, 1 x 16 tubes: 5 seconds 20 $\mu$ L /well, 384 wells, 1 x 16 tubes: 14 seconds 3 $\mu$ L /well 1536 wells, 2 x 32 tubes: 7 seconds  |  |
| Volume range                  | $3$ - 3,000 $\mu L$ /well, selectable in 1 $\mu L$ increments Minimum Prime Volume: 12 mL  |  |
| Flow rates                    | User programmable rates from high to low   |  |
| Dispense accuracy             | ±1 μL at 5 μL and 20 μL; ±1% at 100 μL   |  |
| Dispense precision            | $\leq\!5\%$ CV at 5 $\mu L;$ $\leq\!2.5\%$ CV at 20 $\mu L;$ $\leq\!1\%$ CV at 100 $\mu L$   |  |
| Supply bottle                 | 1 L or 2 L   |  |
| Sterilization                 | Chemical, autoclavable opton   |  |
| Dispensing - RAD              |  |  |
| Labware types                 | Single tip: 6-, 12-, 24-, 48-, 96-, 384-well plates; low profile standard height and deep well formats; PCR trays and microtube 8-to-1 tip: 6-, 12-, and 24-well plates  |  |
| Manifold types                | RAD single, with plastic or steel tip with 1, 5 or 10 tubing, 7° angle RAD 8-to-1 plastic tip, with 5 $\mu$ L tubing, angled bulk dispense chut  |  |
| Volume range                  | 500 nL - 30,000 μL   |  |
| Minimum prime<br>volume       | 1 μL cass, 18": 90 μL ; 1 μL cass, 30": 150 μL<br>5 μL cass, 18": 320 μL; 1 μL cass, 30": 530 μL<br>10 μL cass, 18": 555 μL; 10 μL cass, 30": 920 μL   |  |
| Dispense speed<br>(high flow) | 1 μL cass, 1 μL/well: 19s (96 wells) 55s (384 wells)<br>1 μL cass, 10 μL/well: 33s (96 wells), 112s (384 wells)<br>5 μL cass, 5 μL/well: 19s (96 wells), 58s (384 wells)<br>5 μL cass, 100 μL/well: 76s (96 wells), 286s (384 wells)<br>10 μL cass, 10 μL/well: 21s (96 wells), 66s (384 wells)<br>10 μL cass, 100 μL/well: 70s (96 wells), 259s (384 wells)   |  |
| Dispense<br>performance       | 1 μL cass (med), 0.5 μL/well: Precision 10% CV 1 μL cass (med), 1 μL/well: Accuracy ± 10%, Precision 10% CV 1 μL cass (med), ≥2 μL/well: Accuracy ± 5%, Precision 5% CV 5 μL cass (high), 5 μL/well: Accuracy ± 4%, Precision 5% CV 5 μL cass (high), ≥10 μL/well: Accuracy ± 2%, Precision 2.5% CV 10 μL cass (high), ≥10 μL/well: Accuracy ± 4%, Precision 4% CV 10 μL cass (high), ≥20 μL/well: Accuracy ± 2%, Precision 2 %CV 8-to 1 cass (high), ≥10 μL/well: Precision 2.5%CV 8-to 1 cass (high), ≥80 μL/well: Accuracy ± 2% |  |
| Physical Characte             | ristics  |  |
| Power                         | 100 - 240 Volts AC. 50/60 Hz   |  |
| Dimensions                    | Base instrument: 17.19"W x 11.75" D x 8" H (43.51 x 29.21 x 20.32 cm)  |  |
| Weight                        | Base instrument: 19.5 lbs (8.8 Kg)   |  |
| Regulatory                    |  |  |
| Regulatory                    | CE and TUV marked. ROHS compliant. In Vitro Diagnostic use models are available.   |  |

## MicroFill™ Dispenser

With its microprocessorcontrolled syringe drive technology, the MicroFill™ Microplate Dispenser provides outstanding accuracy and precision while dispensing into 24-, 96- and 384-well plates.

### Low Maintenance Design

The MicroFill is an economical, compact and reliable alternative to conventional microplate dispensers. Its microprocessor-controlled syringe pump provides optimal dispense performance without time-consuming recalibration, cassette replacement and maintenance. Syringes are ideal for higher volume filling, with noteworthy speed improvements compared to other dispense technologies.

### **Guaranteed Sterility**

The entire fluid path is autoclavable for applications requiring sterility. The MicroFill's pump, tubing, dispense manifold and supply bottle are quickly changed for no reagent carryover. User-controlled dispense flow rates allow low- to high-velocity dispensing for both biochemical and cell-based assays. Lowprofile, standard and deep well microplates are all accommodated with a broad volume range from 5 µL to 6 mL.

### **Unattended Operation**

For increased throughput, the MicroFill can be integrated with BioTek's BioStack™ Microplate Stacker or interfaced to third party automated systems with its available interface software. MicroFill drivers are available from most of today's leading system providers.

- Primary and secondary screening assays
- Compound storage
- ▶ Genomics and proteomics research
- Cell-based assays
- ► ELISAs



| Sonoral                     |   |  |
|-----------------------------|---|--|
| General                     |   |  |
| Microplate types            | 24-, 96- and 384-well<br>Low profile, standard and deep well formats  |  |
| Other labware supported     | PCR tubes, microtubes   |  |
| Onboard software            | Create, edit or run multiple protocols  |  |
| Software (pc control)       | Interface software (optional) for robotic system integration  |  |
| Shaking                     | Programmable up to 60 minutes<br>Slow, medium, fast or variable   |  |
| Soaking                     | Programmable up to 60 minutes   |  |
| Automation                  | BioStack™ and 3rd party automation compatible   |  |
| Dispensing - syringe pump ( | (multi-channel)   |  |
| Manifold types              | 924-well dispensing: One 8-tube (1x8) manifold - 316 stainless steel tubes 96-well dispensing: One 8-tube (1x8) manifold - 316 stainless steel tubes 96-/384-well dispensing: |  |
| Dispense speed              | 96 wells, 10 $\mu$ L/well, 1 $\times$ 16: 4 sec 384 wells, 5 $\mu$ L/well, 1 $\times$ 16: 7 sec   |  |
| Volume range                | 5 - 6000 µL/well (manifold dependent)<br>Minimum Prime Volume: 10 mL  |  |
| Flow rates                  | User programmable rates from high to low  |  |
| Dispense accuracy           | $\pm 1~\mu L$ at 5 $\mu L$ and 20 $\mu L$ $\pm 1\%$ at 100 $\mu L$  |  |
| Dispense precision          | ≤5% CV at 5 μL<br>≤2.5% CV at 20 μL<br>≤1% CV at 100 μL   |  |
| Supply bottle               | 1L  |  |
| Sterilization               | Autoclave, chemical   |  |
| Physical Characteristics    |   |  |
| Power                       | 100 - 240 Volts AC. 50/60 Hz  |  |
| Dimensions                  | 18" D x 15" W x 7"H (46 cm x 38 cm x 18 cm  |  |
| Weight                      | 20 lbs (8.9 kg)   |  |
| Regulatory                  |   |  |
| Regulatory                  | CE and TUV marked. ROHS compliant. In Vitro Diagnostic use models are available.  |  |
|                             |   |  |

# Dispenser Comparison Chart

| Which Dispenser is right for you?             | EL406™                     |
|---|----------------------------|
| Key Features                                  |                            |
| ELISA   | •                          |
| Cell-based assays                             | •                          |
| Number of reagents                            | 1 to 3                     |
| Dispense technology                           | Peristaltic and/or Syringe |
| Fully modular and upgradable                  | •                          |
| Performance Specifications                    |                            |
| Dispensing speed                              |                            |
| Peristaltic pump (8-tip, 1x8)                 |                            |
| 96-well, 10 μL/well; 384-well, 5 μL/well      | 3 sec; 6 sec               |
| Syringe pump (16-tube, 1x16)                  |                            |
| 96-well, 20 μL/well; 384-well, 20 μL/well     | 5.25 sec; 14 sec           |
| Dispense accuracy - typical at 5 µL           |                            |
| Peristaltic pump                              | ±2%                        |
| Syringe pump                                  | ±1 μL                      |
| Dispense precision - typical at 5 µL          |                            |
| Peristaltic pump                              | ≤2.5% CV                   |
| Syringe pump                                  | ≤5% CV                     |
| General Specifications                        |                            |
| Microplate types                              | 96, 384 and 1536           |
| Low profile and standard height               | •                          |
| Deep well                                     |                            |
| Strips and full plates                        | •                          |
| Cassette/manifold                             |                            |
| RAD cassettes                                 |                            |
| Peristaltic pump                              | 8-tip (1x8)                |
| Syringe pump, 6-well to 384-well dispensing   |                            |
| Syringe pump, 96-well to 1536-well dispensing | 8-tube (1x8)               |
| Automation ready/BioStack™ compatible         | •                          |
| Variable flow rates                           | •                          |
| Volume range                                  | 1 - 3,000 μL/well          |
| Microplate shaking                            | •                          |
| Autoclavable fluid path                       | •                          |
| Onboard software included                     | •                          |
| Liquid Handling Control™ Software compatible  | •                          |

| MultiFlo <sup>™</sup> FX   | MicroFill™        |
|----------------------------|-------------------|
| Marin 10                   |                   |
| •                          | •                 |
| •                          | •                 |
| 1 to 4                     | 1                 |
| Peristaltic and/or Syringe | Syringe           |
| •                          |                   |
|                            |                   |
|                            |                   |
|                            |                   |
| 3 sec; 6 sec               |                   |
|                            |                   |
| 5.25 sec; 14 sec           | 4 sec; 7 sec      |
|                            |                   |
| ±2%                        |                   |
| ±1 μL                      | ±1 μL             |
|                            |                   |
| ≤2.5% CV                   |                   |
| ≤5% CV                     | ≤5% CV            |
|                            |                   |
| 6 to 1536                  | 24, 96 and 384    |
| •                          | •                 |
| •                          | •                 |
| •                          | •                 |
|                            |                   |
| •                          |                   |
| 8-tip (1x8)                |                   |
| •                          | 8-tube (1x8)      |
|                            | 8-tube (1x8)      |
| •                          | •                 |
| •                          |                   |
| 1 - 30,000 μL/well         | 5 - 6,000 μL/well |
| •                          |                   |
| •                          | •                 |
| •                          | •                 |
| •                          |                   |
|                            |                   |

# Liquid Handling Control™ Software

Liquid Handling Control™ (LHC™) Software allows MultiFlo™ FX Dispenser, EL406™ Washer Dispenser and 405™ Touch and 405 LS Washer users the convenience of programming important assay-specific protocol requirements in a Windows® environment.

#### **Expanded Versatility**

LHC Software is a powerful yet flexible interface for use with BioTek's microplate dispensers and washers. Any programming sequence possible onboard the liquid handler may be duplicated from the computer with LHC Software. LHC also allows a virtually unlimited number of methods to be linked together for the most complex liquid handling routines. From a washer's first prime routine, multiple microplate processes over time, ultrasonic cleaning to dissolve protein or salt crystal build-up to a final system rinse, LHC Software enables unattended operation.

### 21 CFR Part 11 Compliance

To meet the demands of the GxP laboratory, LHC Secure offers features to ensure compliance to 21 CFR Part 11. Flexible multi-user permission levels and electronic protocol and system audit trail signing are all available whenever additional security is required.

#### **Custom Maintenance Reminders**

To facilitate maintenance and keep a washer or dispenser in peak condition, factory recommended maintenance procedure reminders can be preset and customized appropriately for a busy laboratory's usage and throughput requirements. LHC also supports BioStack™ Microplate Stacker and BioSpa™ Automated Incubator integrations.


### Safe Record Keeping

Protocol parameters may be quickly printed for safe record keeping. Alternatively, onboard instrument protocols may be uploaded and backed up on a laboratory's network. Satellite research labs working on joint projects can be certain their wash parameters are identical for experimental integrity.

### SiLA Compliant Drivers

For automated systems that require SiLA compliant integration, LHC SiLA is available.





# Washer & Dispenser Accessories

BioTek offers a wide range of accessories to help increase productivity, expand your plate washer's and dispenser's capabilities, and maintain the performance of your BioTek microplate liquid handling system. See our web site for a complete listing of available accessories.



### **Peristaltic Pump Dispenser Cassettes**

A wide selection of peristaltic pump cassettes are available with choices in volume ranges, tip materials and bore sizes for use with  $\mathsf{EL406}^{\mathsf{TM}}$  and  $\mathsf{MultiFlo}^{\mathsf{TM}}$  FX.



**BioStack™ Microplate Stacker** 

Automate routine processes with this compact stacker.



BioSpa<sup>™</sup> 8 Automated Incubator

BioSpa 8 is an automated incubator linking BioTek readers or imagers together with washers and dispensers for full workflow automation of up to 8 microplates.



### Dispense/Waste Systems

A dispense/waste system is required on all 405™ Touch, 405 LS and EL406 models. Many selections are available based on throughput, bottle size and vacuum pump requirements.



### Syringe Pump Dispenser Manifolds

A range of MultiFlo FX and MicroFill™ dispense manifolds are available for various microplate types and reagent characteristics.



#### 3-Instrument Rack

For third party robotic system integration, a rack is available for supporting up to three dispensers or other BioTek instrumentation.



### Instrument Qualification

See the Compliance Section on pages 62-63 for details about BioTek's product qualification tools and services.

# Precision™ Microplate Pipetting Systems

The Precision™ is an innovative solution for automated liquid handling. With its ability to perform virtually any routine liquid transfer, Precision replaces tedious manual pipetting.

### **Automate Manual Pipetting**

The Precision can be customized with a range of options perfect for medium throughput labs looking to automate their everyday pipetting with walk-away confidence. BioTek's proprietary pipette technology and unique tip sealing allow most standard tips to be used for transfers in common sample tubes and 6- to 384-well microplate formats.

# Open Deck Layout and Flexible Software

A user-configurable, multistation deck allows for flexible experimental design; microplates, tips and other labware may be placed in nearly any location for optimal efficiency. Available Precision Power™ Software offers complete Precision control with intuitive protocol creation, expanding the instrument's dynamic capabilities with a graphical program simulator and sample ID tracking.

# Space Saving, Compact Footprint

Its small footprint and wellorganized design make the Precision ideally suited for installation inside standard size biological safety cabinets and chemical fume hoods. The Precision XS model delivers outstanding liquid handling performance with four liquid transfer tools on a single platform. All four may be intermixed throughout a fully automated protocol - single and multichannel pipetting along with single- and multi-channel bulk reagent dispensing.

- ▶ Sample transfers from tube to microplate
- Serial dilutions
- Mixing
- ▶ Plate replication mother/daughter transfers
- Reagent addition
- Hit picking
- ELISA automation
- Secondary screening assays
- ▶ Compound profiling
- Cell-based assays



| General                    | General  |  |  |
|----------------------------|--|--|--|
| Microplate types           | Precision: Occided 284 well microplates,   |  |  |
|                            | 96- and 384-well microplates   |  |  |
| Other labware<br>supported | Test tubes <100 mm   |  |  |
| Onboard software           | Precision:Create, edit or run multiple protocols<br>Precision XS: Computer control only  |  |  |
| Software<br>(pc control)   | PrecisionPower Software  |  |  |
| Automation                 | BioStack™ and 3rd party automation compatible (except 1 x 12 pipetting)  |  |  |
| Platform                   | Precision XS: 6 stations<br>Precision: 1x8, 6 stations; 1x12, 4 stations   |  |  |
| Dispensing - Syringe Pump  |  |  |  |
| Manifold types             | Precision XS: 96-/384-well dispensing: One 8-tube (1x8) manifold - 316 SS tubes 6- to 384-well dispensing: One single-channel probe  Precision: One 8- (1x8) and/or 12-tube (1x12) manifold - 316 SS tubes |  |  |
| Dispense speed             | Precision XS: 100 µL/well, 96 wells, 1x8: 14 seconds 100 µL/well, 96 wells, single-channel: 4 minutes  Precision: 100 µL/well, 96 wells, 1x8: 14 seconds   |  |  |
| Volume range               | Precision XS: 1x8: 10 $\mu$ L - 10 mL Single-channel: 5 $\mu$ L - 10 mL Precision: 10 $\mu$ L - 10 mL  |  |  |
| Dispense accuracy          | ±1% at 100 μL  |  |  |
| Dispense precision         | ≤1.5% CV at 100 μL   |  |  |
| Supply bottle              | 2 L; 2 L and 125 mL (Precision XS)   |  |  |
| Sterilization              | Autoclave, chemical  |  |  |
|                            |  |  |  |


| Pipetting                |  |  |
|--------------------------|--|--|
| Manifold types           | Precision XS: 1x8, single channel<br>Precision: 1x8, 1x12 (configuration dependent)  |  |
| Pipetting speed          | Precision XS:<br>1x8: 100 µL/well, 96 wells, tip change: 3 min<br>Single-channel: 100 µL/well, 96 wells, tip change: 22 min<br>Precision:  |  |
|                          | 1x8: 100 μL/well, 96 wells, tip change: 3 min  |  |
| Volume range             | Precision XS:<br>1x8: 5 - 120 μL<br>Single channel: 5 - 200 μL   |  |
|                          | Precision: 5 - 120 μL  |  |
| Fluid delivery           | Air displacement syringe drives  |  |
| Dispense accuracy        | ±1% at 100 μL  |  |
| Dispense precision       | ≤1.5% CV at 100 µL   |  |
| Pipette tips             | BioTek and other commercially available tips   |  |
| Physical Characteristics |  |  |
| Power                    | 100 - 240 Volts AC. 50/60 Hz   |  |
| Dimensions               | Precision XS Instrument only: $25^{\circ}$ W × $16^{\circ}$ D × $20^{\circ}$ H ( $640$ × $410$ × $510$ mm ) Instrument with optional cabinet: $30^{\circ}$ W × $20^{\circ}$ D × $22^{\circ}$ H ( $760$ × $510$ × $560$ mm)  Precision Instrument only: $21^{\circ}$ W × $15^{\circ}$ D × $16^{\circ}$ H ( $525$ × $374$ × $400$ mm) Instrument with optional cabinet: $23^{\circ}$ W × $17.5^{\circ}$ D × $17.5^{\circ}$ H ( $575$ × $435$ × $435$ mm) |  |
| Weight                   | Precision XS: Instrument only: 40 lbs (18 kg) Instrument with optional cabinet: 64 lbs (29 kg) Precision: Instrument only: 28 lbs (12.7 kg) Instrument with optional cabinet: 38 lbs (17.2 kg)   |  |
| Regulatory               |  |  |
| Regulatory               | CE and TUV marked. ROHS compliant. In Vitro Diagnostic use models are available.   |  |
|                          | I  |  |





"This stacker is very simple to program and operate through a PC based software. It is a very reliable machine. Its open format allows for easy cleaning and maintenance and the release mechanism is up front and easily accessible for rapid removal and installation of the stacks."

(BioStack Microplate Stacker)



Automation products provide speed, flexibility, consistent results and unattended operation when configured with BioTek's line of microplate imagers, readers, washers and dispensers. The result is a scalable, cost-effective system that can adapt to your changing requirements. The latest innovation in BioTek's line of automation products, the BioSpa™ 8 Automated Incubator, is designed to automate cell-based assay workflows in a controlled and monitored environment.

### BioSpa<sup>™</sup> 8 Automated Incubator

BioTek's new BioSpa™ 8 automates incubated assay workflows by moving and storing microplates containing live cells or temperature sensitive reagents. More versatile than a benchtop incubator, BioSpa 8 manages up to 8 microplates, flasks, or cell culture dishes in a CO<sub>2</sub>/O<sub>2</sub>, temperature and humidity monitored environment. Integrated with BioTek's washers, dispensers, imaging and detection systems, BioSpa 8 manages the entire process from sample preparation to detection or imaging in one compact system.

# Environment control and monitoring leads to cell assay success

BioSpa 8 offers incubation to 45 °C,  $CO_2/O_2$  control and monitoring, plus humidity monitoring – everything a successful live cell assay needs.

### Biosafety cabinet compatible

BioSpa 8 is designed to help protect against contamination, with a HEPA filter for incoming air and an interior that is easily cleaned and decontaminated. For the ultimate protection against potentional contamination, BioSpa 8 is compact – it fits within a biosafety cabinet along with integrated washer, dispenser, imager or plate reader.

# Full workflow automation integrates sample prep

BioSpa 8 automates processes that commonly burden many labs working with live cells; inconvenient culture maintenance requirements, contamination hazards and handling multiple instruments required for both sample plate preparation and downstream processing. BioSpa 8 handles from 1 to 8 plates, moving them between the integrated washer or

dispenser and imaging system or multi-mode reader for complete, unattended process automation.

# Continuous recording and monitoring with notifications

BioSpa 8 continuously monitors and records important workflow parameters, and can automatically send text or email notifications. BioSpa 8 provides confidence and control for unattended automation.

# Simple integration for rapid implementation

BioSpa 8 is compatible with several BioTek imaging and multimode readers, plate washers, dispensers and combination systems. The simple integration doesn't require specialized tables or other hardware or software, and BioSpa 8 is compact enough to be used in a biosafety cabinet for critical live cell assays.

### Typical Applications:

Automated sample preparation for cell based assays


Drug Absorption

Automated sample preparation for cell based assays

- Cell Culture QC
- Cell Proliferation
- Apoptosis
- Cytotoxicity
- 3D Cell Culture
  - Tumor Invasion
  - Signal Transduction

  - Phenotypic Assays
- ▶ Cell Migration and Invasion
- ► Fluorescent Protein Detection
- RNA Expression





| General                             |  |  |
|-------------------------------------|--|--|
| Microplate types                    | 6- to 1536-well standard height microplates, with or without lids. Plate height range: 7.6 mm to 25.4 mm   |  |
| Other labware supported             | Petri and cell culture dishes (35 mm and 60 mm), T25 flasks  |  |
| Plate capacity                      | Up to 8 microplates  |  |
| Air filter                          | User-replaceable HEPA filter   |  |
| Decontamination:                    | Easy interior access for cleaning and decontamination.   |  |
| Plate handler                       | Robotic arm moves plate to and from BioSpa 8 and connected instrument; handles de-lidding and re-lidding plates.   |  |
| Dimensions                          | 27.2" W x 20.1" D x 18.9" H<br>69 cm W X 51 cm D x 48 cm H   |  |
| Compatible BioTek instruments       | Cytation 5, Cytation 3, Synergy Neo2, Synergy H1, Epoch 2<br>EL406, 405 TS, 405 LS, MultiFlo FX, MultiFlo  |  |
| Interfacing capacity                | 1 or 2 devices: Reader/imager only, washer/dispenser only, or both   |  |
| Temperature Control                 |  |  |
| Range                               | To 45 °C   |  |
| Control resolution                  | 0.1 °C   |  |
| Uniformity                          | ± 0.5 °C at 37 °C  |  |
| CO <sub>2</sub> Control             |  |  |
| Range                               | 1 – 19%  |  |
| Control resolution                  | ±0.1%  |  |
| Stability                           | ±0.2 at 1% O <sub>2</sub>  |  |
| O <sub>2</sub> Control              |  |  |
| Range                               | 1 – 19%  |  |
| Control resolution                  | ±0.1%  |  |
| Stability                           | ±0.2 at 1% O <sub>2</sub>  |  |
| Humidity                            |  |  |
| rH                                  | 80 to 95% (lidded plates and 5% CO <sub>2</sub> )  |  |
| Source                              | Removable water pan  |  |
| Water level sensor                  | Low water level alert  |  |
| Software                            |  |  |
| BioSpa Automated Incubator Software | Provides programming interface for BioTek detection and liquid handling devices Allows user notification (text or email) of events and/or errors in the system Provides control, monitoring and logging of: CO/O2 - temperature control Provides humidity level monitoring and logging Allows long-term uninterrupted runs up to 2 weeks |  |
| Regulatory                          |  |  |
| Regulatory                          | CE and TUV marked, RoHS compliant.   |  |

## BioStack™ Microplate Stacker

BioStack™ is a compact and versatile microplate stacker compatible with BioTek's microplate washers, dispensers, pipetting, detection and imaging systems. BioStack is easy to use and provides walkaway automation for routine processes, including processes requiring plate de-lidding and re-lidding.

### **Ultra Fast Transfer Speeds**

BioStack offers the fastest plate transfer time, taking less than 10 seconds to remove and replace plates on the instrument carrier. BioStack is well-suited for high throughput plate stacking requirements with BioTek readers, washers and dispensers.

### Plate De-lidding

Many cell-based microplate processes require lidded plates during incubation and to protect sterility. Typically, automation of these processes meant purchasing an expensive microplate handler to de-lid the plates for measurement or liquid handling operations. BioStack now offers an affordable option for plate de-lidding in the BioStack 4 model to interface with BioTek's detection and liquid handling instruments.

### Multiple Microplate Geometry Compatible

BioStack is compatible with standard 96- and 384-well plates, low volume 384-well plates and 1536-well plates. The BioStack 4 adds 24- and 48-well plates to its menu of compatible microplate labware, providing higher throughput in a walk-away system for a variety of microplate geometries. An available barcode scanner provides additional automation for high-throughput plate processing.

Plate IDs are read and sent to the plate data file in Gen5<sup>™</sup> or LHC<sup>™</sup> Secure software for storage or export.

#### 10-, 30- or 50-Microplate Stacks

Choose between 10-, 30- or 50-plate stacks to best suit your throughput requirements. Low volume, half-height plates are also compatible, with up to 75-plates capacity in the 50-plate stack.

### Compact, Rugged Design

BioStack allows worry-free operation, even under the heaviest usage. The motors, mechanical assemblies and software are all designed for long term, continuous use and maintenance-free use. The rotational gripper and very small footprint allows for integration position versatility and for optimal fit within a biosafety enclosure or for space-savings on the benchtop.


- Cell-based assays
- ▶ ELISAs
- Primary screening assays
- Colorimetric, fluorometric and luminescent assays





| BioStack 4 | BioStack 3 • BioStack Neo • BioStack |
|------------|--------------------------------------|
|            |                                      |

| General                                |  |   |  |
|--|--|---|--|
| Microplate types                       | ANSI/SLAS Standard and low profile 96-, 384- and 1536-well plates 24- and 48-well plates (model dependent) Maximum plate height 23.2 mm  | ANSI/SLAS Standard and low profile 96-, 384- and 1536-well plates Maximum plate height 14.6 mm  |  |
| Lidded plate handling                  | De-lidding capability: (lids always removed during processing) 96-,384-, and 1536-well plates. Maximum height, including lids:16.9 mm  Nunc plates: (lids can remain on plate during process, or can be removed) 6-, 12-, 24-, 48- well plates.  Maximum height, including lids: 23.2 mm | n/a   |  |
| Microplate capacity                    | 10 and 30 plate stacks are removable and interchangeable (50-plate stacks may be used with non-lidded plates only) 96-/384-well plates: Up to 30 plates (with lids) 1536-well plates: Up to 75 plates  | 10-, 30- and 50-plate stacks are removable<br>and interchangeable<br>96-/384-well plates: Up to 50 plates<br>1536-well plates: Up to 75 plates                                |  |
| Barcode scanner<br>(option)            | Landscape or portrait orientation, Code 39, Codabar, UPC/EAN, Code 128 compatible  | Landscape or portrait orientation, Code 39, Codabar, UPC/EAN, Cod<br>128 compatible<br>1D and 2D barcodes (BioStackNeo)   |  |
| Processing speed (plate exchange time) | <20 seconds (with de-lidding)<br><12 seconds (without lids)  | <10 seconds: BioStack 3, BioStack Neo<br><33 seconds: BioStack  |  |
| Direct control                         | Washers and dispensers with keypad interface can directly control BioStack   | Washers and dispensers with keypad interface can directly control BioStack  |  |
| PC software                            | LHC2 for liquid handling instruments (optional)<br>Gen5 2.0 for readers  | LHC2 for liquid handling instruments (optional)<br>Gen5 2.0 for readers   |  |
| Physical Characteristi                 | ics  |   |  |
| Power                                  | 100 - 240 Volts AC. 50/60 Hz.  | 100 - 240 Volts AC. 50/60 Hz.   |  |
| Weight                                 | <25 lbs (11.3 kg)  | <25 lbs (11.3 kg)   |  |
| Dimensions                             | 8.3" W x 22" D (21 x 56 cm)  Overall height will vary depending on connected instruments and stacks used   | BioStack and BioStack Neo 7.4" W x 20.7" D (18.8 x 52.6 cm)  BioStack 7" W x 18.5" D (18 x 47 cm) Overall height will vary depending on connected instruments and stacks used |  |
| Regulatory                             | '  |   |  |
| Regulatory                             | CE and TUV marked. ROHS Compliant.   |   |  |



"BioTek is great. The BioTek engineer is very professional, helpful, knowledgeable, and flexible."



Our teams are committed to providing the service and support you need to sustain the optimal performance of your BioTek products. BioTek Service Engineers provide personal support for instrumentation, software, parts and applications at our Global Technical Support Center. BioTek Scientists, Engineers, Technicians and Sales Representatives provide valuable assistance to laboratories worldwide.

### Compliance

As an ISO certified manufacturer, BioTek understands the importance of standardized product qualification procedures and traceability and provides a number of tools and services designed to streamline the process and minimize the resources required to perform such testing.



# 21 CFR Part 11 Compliance Products

Gen5™ Secure Data Analysis
Software for microplate readers and
LHC™ Secure protocol definition
and control software for BioTek's
microplate washers and dispensers
are uniquely designed to help
ensure compliance to 21 CFR
Part 11. Both software programs
offer important security features,
including:

- ► Electronic signature of data and protocol files
- Secure data storage
- Multiple and definable user permission levels
- ▶ Data and protocol audit trails
- ▶ Protected functions

#### **IVD Compliance**

Many BioTek microplate instruments are labeled for In Vitro Diagnostic use, identified by the IVD logo.
Other products may have IVD



Compliant models available. Contact

CustomerCare@biotek.com for more information.

### **Product Qualification**

### **Software Validation**

A Validation Package is available for Gen5 Software to allow testing and validation of key functions within Gen5 and Gen5 Secure. Included in the easy-to-use package are:

- ► Test Plans
- ▶ Results Checklists
- ► Data sets



### Instrument IQ/OQ/PQ Packages

BioTek offers a complete menu of Product Qualification Packages for all of our microplate instruments. All product Qualification Packages are fully validated to assure that the procedures and associated data/spreadsheets supplied in the package meet regulatory requirements. Within each package, you'll find detailed:

- ► Product Specifications
- Qualification interval guidelines
- ► IQ/OQ/PQ test plans and procedures
- ▶ Data sets (where applicable)
- Qualification checklists and log sheets for complete documentation

### RoHS2 Directive 2011/65/EU

BioTek is committed to helping protect the environment in all of our customers' countries. BioTek products that meet the RoHS directive are indicated in the Regulatory section of the product specifications in this catalog.





### **Test Plates**

The use of standardized plates to supplement the verification of an instrument's performance is a time-and resource-saver in most laboratory environments. BioTek offers several test plates to facilitate the test procedures found in our microplate reader IQ/OQ/PQ packages, and can be automated through the Gen5 software.

#### **Absorbance Test Plates**

For use with the ELx800, ELx808, Epoch 2, Epoch, PowerWave HT, Synergy and Cytation multi-mode reader absorbance modes. Ensure GxP compliance by checking instrument performance against specifications for:

- Accuracy
- ► Repeatability
- ▶ Linearity
- Wavelength accuracy (for monochromator-based systems)

#### Fluorescence Test Plates

Ideal for quick checks of the fluorescence intensity detection system between more thorough instrument qualification. The Fluorescence Test Plate aids in maintaining GxP compliance by automatically checking a series of critical performance parameters, including:

- ► Alignment
- ► Cross talk
- ► Signal-to-noise ratio
- ▶ Linearity
- ▶ Precision

### **Luminescence Test Plates**

This NIST-traceable Luminescence Test Plate is used with the applicable Product Qualification Package or updated User's Manual. Features include:

- ➤ NIST-traceability certificate guarantees a controlled light output from the test plate
- ➤ Simple design, easy to use: just turn the plate on, and read the ultra-stable, low light level LEDs

Test Plate Recertification Programs are available. Contact BioTek Service for details. www.biotek.com/contact

### **Applications Support**

At BioTek, our customers' applications come first. Whether it's our existing or new products in development, we design in capabilities to enable your most important applications. To facilitate this process, BioTek has an on-site Applications Lab with a team of seasoned scientists continuously working on cutting-edge scientific applications and partnering with the best-known reagent and consumable vendors. Below are some examples of hot application areas where we have demonstrated the utility of our products.



#### 3D Cell Culture

- Signal transduction in a collagenbased scaffold
- Long term toxicity in liver microtissues
- Methods development for spheroid formation in hanging drop plates
- Automation of 3D cell culture work flows
- Cell invasion assays using spheroids formed in ULA microplates

#### Phenotypic Assays and Screening

- Hypoxia assays in keratinocytes and spheroids
- Oxidative stress assays for ROS production
- Cell cycle using nuclear stains and sensors

- Mitochondrial oxidative stress and apoptosis
- Autophagy and lysosomal disorder assays

### Live Cell Assays

- RNA quantification using fluorescent nanoprobes
- Multiplexed second messenger assays (Ca<sup>2+</sup>, cAMP, diacylglycerol, PIP<sub>2</sub>) using genetically engineered probes
- ► Cell invasion and migration in FluoroBlok microplates
- Automation and analysis of drug absorption in Caco-2 and MDCK cells
- Characterization of multi-drug resistance transporters

### Biologics/Biosimilars

- ADCC assays using nonradiometric detection and freshly isolated NK cells, cryopreserved NK cells, an NK cell line and a bioluminescent cell reporter
- Immunogenicity assays with AlphaLISA beads
- Bridging assay comparisons between solution ELISA and AlphaLISA

- Cell-based assays for blocking antibodies
- Aggregation assay using a molecular rotator-type fluorescent probe

### **Food Safety and Quality**

- Food freshness assay based on ATP depletion
- Determination of E coli and other pathogens in lettuce wash
- Antioxidant potential using an ORAC Assay
- Semi-automated ELISA assay for melamine in milk
- ▶ Analysis of histamine in wine

### **Biofuel Research**

- Determination of algal cell lipids using Nile Red
- Monitoring algal growth using their intrinsic properties
- Enzymatic digestion of polysaccharides
- Identification of biofuel producing bacteria through temperature resistance
- Monitoring enzymatic glucose production from cellulosic feedstock

More than 2,600 BioTek application notes, white papers, poster presentations, citations and sample files are available at:

www.biotek.com/resources

### Global Service & Support

Extend the life of your BioTek instrument, and protect your research results, with BioTek's service professionals. Our service experts in the field and at our regional service centers receive extensive, ongoing training at our headquarters to stay abreast of the latest products, and service techniques. Our products and services are compliant with FDA, GLP and ISO requirements. With all of this information at hand, our service experts help you to maintain precise results over the life of your BioTek instrument while providing an experience that is superior to our competition. For any service or support need, contact us at TAC@biotek.com or (888) 451-5171.



#### Field Service

Our team is ready to visit your laboratory and provide:

- Installation, Training and Installation Qualification
- Operational Qualification
- Preventive Maintenance
- Instrument Upgrades and Software Upgrades
- Repairs

### **Regional Service Centers**

BioTek Service Centers are located across the globe ready to service your BioTek products:

- ► Test Plate Certification
- Preventive Maintenance
- Instrument Upgrades
- ▶ Dispense Cassette Refurbishment

### Technical Assistance Center (TAC)

BioTek's TAC is staffed with skilled scientists and engineers available to provide technical assistance for instrumentation, software and applications.

#### **Customer Resource Center (CRC)**

BioTek's Customer Resource Center gives customers access to information about their specific BioTek microplate instrumentation and software. This web site makes it easy for customers to acquire relevant and necessary information about their products.

#### Customers can:

- Track orders
- Maintain equipment inventory
- Access warranty information
- Download technical information, user manuals and software updates
- Request service and technical support



BioTek Preventive Maintenance Service includes a certificate of calibration for every instrument.

# Third party customer satisfaction survey excerpts:

"I was discussing with a colleague after the help BioTek provided how competent and willing to help BioTek has always been. Because you guys are willing to just sort out problems, even with old equipment, you have a loyal customer base. I have and will continue to buy BioTek products wherever I work. So, thanks, we look forward to continued positive interactions with BioTek."


"BioTek provided the best customer service. Service was able to listen to my problem and resolve it very quickly. I am very pleased by the level of service I received. I actually enjoy calling BioTek"

"Working with BioTek is really Awesome. I get my answer quickly, A high level of follow up. Always good results. Just keep it like this."

details on BioTek's service and support are available at www.biotek.com

Access to BioTek's Customer

Resource Center and more



BioTek is a family-run organization founded in 1968. In 1981 BioTek entered the microplate instrument arena and introduced its first microplate reader. Since then, BioTek has emerged as a global leader in microplate-based solutions that increase the productivity for customers engaged in healthcare, pharmaceutical, agricultural and research applications. Today, BioTek is completely focused on microplate instrumentation, automation and software. Our products continue to be designed and manufactured in Vermont, U.S.A.

1968 - 79 < BioTek develops, manufactures and sells biomedical testing equipment

**1969** BioTek moves to its first facility located on East Spring Street, Winooski, VT

1980 < BioTek's first laboratory instrument, the EL307 EIA is launched

> BioTek moves to Burlington, VT

1983 < BioTek awarded U.S. patent for EL307's microwell position indicator

▶ 1982 EL307B microplate reader introduced by



1987 < BioTek introduces the EL320 Stacking Automated Reader for automatic reading of up to 25 plates

▶ 1988 BioTek introduces the EL311 and EL312 Microplate Readers, EL403 Automated Microplate Washer, the EL301 Manual Strip Reader for field use, and Kineticalc PC Software



1968 - 1979 1980 - 1985 1986 - 1990

> 1984 BioTek's first Automated Microplate Reader, the EL310, the EL308 Microplate Reader and the first Microplate Washer, the EL402, are introduced

▶ 1985 BioTek introduces the EL309 Microplate Reader, and the EL401 Microplate Strip Washer

BioTek moves to current location in Winooski, VT

**EL340 Biokinetics** reader, EL944 Turbo Software for EL311

1989 <

and EL312 clinical data reduction and the ELP-35 Automated Strip Washer are introduced by BioTek



#### 1991 < Kineticalc II software is introduced by

BioTek

BioTek begins development of Omni System for automated high throughput infectious disease diagnostics and screening



### ▶ 1992

Ceres 900 launched by BioTek

BioTek enters into an agreement with Immucor, Inc, to develop the first automated blood typing and crossmatching analyzer, the ABS2000



#### 1996 <

BioTek receives ISO9001 certification



BioTek introduces ELx808

### 1998 < μQuant and PowerWaveX Microplate Spectrophotometers and KCJunior software are



introduced by

BioTek

### ▶ 1997

BioTek introduces the PowerWave, our first Microplate Spectrophotometer, the FL600 Multi-Detection Reader, the ELx50 Microplate Strip Washer and KC4 Data Reduction Software



#### 2001 <

Powerwave HT, MicroFill are launched by BioTek



#### ▶ 2002

BioTek sells Biomedica division to focus solely on Laboratory Micropla

BioTek opens European Coordination Center in Germany

BioTek launches Synerg HT, PowerWave XS, BioStack



#### 1996 - 2000 2001 - 2005

BioTek launches ELs1000 automated ELISA system

### 1995 4

ELx800 Microplate Reader, EL404 Microplate Washer, and KC3 for Windows data analysis software are introduced by BioTek



### ▶ 1994

BioTek introduces the FL500 Fluorescence Microplate Reader and ELP-40 Microplate Strip Washer



### 1999 <

BioTek launches FLx800, ELx405, ELx405M



BioTek awarded U.S. patent for the ELx405 "Dual-Action" manifold

BioTek awarded U.S. patent for the quartz BioCell cuvette for fixed 1 cm vertical photometry

▶ 2000

Precision 2000 Microplate Sample Processor and Precision Power Software introduced by BioTek



#### 2003 <

BioTek enters into an agreement with UVM and Immucor, Inc to develop a practical fluorescence based platelet assay for determining platelet activation

### ▶ 2004

Precision XS, Clarity introduced by BioTek



### 2005 <

BioTek launches its customer newsletter, TekTalk

BioTek finalizes an agreement with Immucor to develop the next generation blood typing and cross matching instrument, the Galileo Echo

#### 2006 < Synergy 2, NanoQuot, Gen5 software are launched by BioTek



▶ 2007

offices

BioTek opens China,

India, Singapore

Synergy 4, Liquid

Handling Control software are

BioTek launches online Customer Resource Center

introduced by BioTek

2008 4 BioTek launches MicroFlo Select, EL406



2011 < BioTek opens South Korea office

BioTek introduces ELx405 Select Deep Well, Eon, Gen5 version 2.0



### ▶ 2012

405 Touch, Synergy Neo, BioStack3, Gas Controller for Synergy H1, 405 LS are launched by BioTek



BioTek opens Shanghai office

Web site offers content in 10 languages

2015 < BioTek launches Synergy Neo2



BioTek introduces BioSpa 8



BioTek expands in the UK

2011 - 2014 2006 - 2010

Epoch, Take3, Synergy Mx are introduced by BioTek



### BioTek is awarded Best Place to Work in Vermont from VT Business Magazine and Business

Champlain Regional Chamber of Commerce BioTek launches multi-

of the Year from the Lake

lingual web site BioTek receives EP Patent for Verify technology and Ultrasonic Advantage

washer technologies

### ▶ 2010 BioTek launches Synergy H1, Synergy H4, Synergy 2 Alpha, MultiFlo



BioTek wins Vermont Deane C. Davis **Outstanding Business** of the Year award

BioTek opens Switzerland office

### 2013 <

BioTek introduces MultiFlo FX, BioStack 4 and 405 Verify



BioTek opens Taiwan and Japan offices

BioTek enters the imaging market with the Cytation5 Cell Imaging Multi-Mode Reader



BioTek launches MultiFlo FX RAD, Epoch 2, Synergy HTX and Cytation 5



BioTek receives US Patent for Verify technology and Ultrasonic Advantage washer technologies



#### World Headquarters, United States

BioTek U.S.

Phone: (802) 655-4740 Toll Free: (888) 451-5171

Service Toll Free: (800) 242-4685 Email: CustomerCare@biotek.com Service Email: TAC@biotek.com

www.biotek.com

#### China

BioTek China-Beijing Phone: +86 (10) 85865569 Email: infochina@biotek.com www.biotekchina.com.cn

BioTek China-Shanghai Phone: +86 (021) 50435800 Email: infochina@biotek.com www.biotekchina.com.cn

### France

BioTek France

Phone: +33 (0) 389206329 Email: info@biotek.fr www.biotek.fr

### Germany

BioTek Germany European Coordination Center Phone: +49 (0) 71369680 Email: info@biotek.de www.biotek.de

#### India

BioTek India

Phone: +91 7718810038 Email: biotek@biotek.in

www.biotek.in

### Japan

BioTek Japan

Phone: +81 (0)3 5812 8109 Email: infojapan@biotek.com

www.biotek.jp

#### South Korea

BioTek South Korea Phone: +82 (0)2 5624740 Email: korea@biotek.com www.biotekinstruments.co.kr

#### **Switzerland**

BioTek Switzerland Phone: +41 (41) 2504060 Email: info@biotek.ch www.biotek.ch

### Taiwan

BioTek Taiwan

Phone: +886 (2) 26277725 Email: infotaiwan@biotek.com

www.biotek.tw

### **United Kingdom**

BioTek UK

Phone: +44 (1793) 575020 Email: infouk@biotek.com

www.biotek.co.uk

For the rest of the world, contact your local distributor

or BioTek U.S. at:

Phone: +1 (802) 655-4740

Email: CustomerCare@biotek.com

<sup>&</sup>quot;Alpha" refers to products/technologies from PerkinElmer, Inc., that carry trademarks or registered trademarks DLR and DLReady logo are trademarks of Promega Corporation.











Think Possible



At BioTek, our philosophy transcends conventional thinking and challenges the old ways. We develop fresh, original solutions by unifying concepts that often appear to be opposed. It means to shape and reshape. To engineer, build, deliver and support products that <u>best</u> serve the marketplace by providing <u>what</u> you need, <u>when</u> you need it.

Think Possible. It's the difference between leading and following.